 **ORGENTEC** Produktkatalog  
product catalogue

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▶ ▶ ▶ The challenges in laboratory diagnostics are global.  
We face them together with our partners.

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### Fueled by Innovation

— With pioneering spirit, nearly inexhaustible supply of ideas, and high technical expertise, ORGENTEC has become firmly established in medical diagnostics. Our employees' talent and skills establish the basis for our success. At our headquarters in Mainz 130 people are busy in research and development, manufacturing and quality control, distribution and sales. Their commitment is the reason why laboratories everywhere rely on the high standards and excellent service provided by ORGENTEC Diagnostika.

— New ideas move us forward and guarantee future innovation. Creative scientists, experienced engineers, and highly educated technicians work with great dedication on the development of new products and the continuous improvement of our established test systems.

— ORGENTEC Diagnostika excels with superior product quality and excellent service. These originate at our in-house production facility at the company headquarters in Germany. We are certain that we only truly control those things we develop and produce under our own management. We also maintain strict control over the complete logistics chain. We are thus able to maintain our own high standards of quality.

— Ideas, corporate courage and steady growth – these form the basis for more than 30 years of success at ORGENTEC Diagnostika. We are contributing to the benefit of patients with cutting-edge research and technological progress.

# Guide to Autoimmune Serology

## Diseases Diagnostic Marker

Diseases	ANA	ENA	C1q	Centromere B	dsDNA	Histone	Jo-1	Nucleosome	Rib-P	RNP/Sm	RNP-70	Scl-70	Sm	SS-A	SS-A-52	SS-A-60	SS-B	ssDNA	DNAse activity	Rheumatoid Factor	ACPA	MCV	CCP	
Rheumatoid arthritis (RA)	Grey	Grey	Grey			Grey								Grey						Red	Red	Red	Red	
Felty's syndrome	Grey	Grey	Grey			Grey													Grey		Red			
Juvenile chronic arthritis	Grey																		Grey		Grey	Grey	Grey	Grey
Sjogren's syndrome	Grey	Grey													Red	Red	Red	Red			Grey			
Systemic lupus erythematosus (SLE)	Red	Grey	Grey		Red	Grey		Grey	Grey	Grey	Grey		Red	Grey	Grey	Grey	Grey	Grey	Grey	Grey				
Lupus nephritis	Red	Grey	Red		Red	Grey		Grey	Grey	Grey	Grey		Red	Grey	Grey	Grey	Grey	Grey	Grey	Grey				
Neonatal lupus erythematosus	Red	Grey									Grey			Red	Red	Red	Red			Grey				
Congenital heart block (CHB)	Red	Grey												Red	Red	Red	Red			Grey				
Subacute cutaneous lupus erythematosus (SCL)	Grey	Grey												Grey	Grey	Grey	Grey		Grey					
Drug-induced lupus	Grey					Grey		Grey											Grey					
Systemic scleroderma	Grey	Grey		Red		Grey		Grey		Grey	Grey	Red												
Raynaud's syndrome	Grey	Grey		Grey						Grey	Grey	Grey									Grey			
CREST syndrome	Grey			Red								Grey												
Polymyositis/dermatomyositis	Grey	Grey					Red							Grey	Grey									
Idiopathic myositis	Grey						Grey			Grey	Grey			Grey	Grey									
Mixed connective tissue disease (MCTD)	Grey									Red	Red													
Undifferentiated connective tissue disease (UCTD)	Red	Grey								Grey	Grey			Grey	Grey	Grey								
Anti-synthetase syndrome	Grey						Red									Grey								
Cryoglobulinemia																					Grey			
Hypocomplementemic urticarial vasculitis syndrome (HUVS)			Grey																					
Antiphospholipid syndrome (APS)																								
Systemic vasculitis																								
Granulomatosis with polyangiitis (GPA)																								
Microscopic polyangiitis (MPA)																								
Goodpasture syndrome																								
Celiac disease																								
Pernicious anaemia																								
Inflammatory bowel diseases																								
Autoimmune hepatitis type I (AIH I)	Red																							
Autoimmune hepatitis type II (AIH II)																								
Primary biliary cirrhosis (PBC)	Grey			Grey																				
Primary sclerosing cholangitis (PSC)	Grey																							
Hashimoto thyroiditis																								

# CONTENTS

## AUTOIMMUNE DISEASE DIAGNOSTICS

### Rheumatology Diagnostics

ANA Detect . . . . .	8
ANAcombi . . . . .	9
ANAscreen . . . . .	10
Anti-alpha-Fodrin . . . . .	11
Anti-C1q . . . . .	12
Anti-CCP <i>hs (high sensitive)</i> ® . . . . .	13
Anti-Centromere B . . . . .	14
Anti-dsDNA . . . . .	15
Anti-Histone . . . . .	16
Anti-Jo-1 . . . . .	17
Anti-MCV® . . . . .	18
Anti-Nucleosome . . . . .	19
Anti-Rib-P . . . . .	20
Anti-RNP/Sm . . . . .	21
Anti-RNP-70 . . . . .	22
Anti-Scl-70 . . . . .	23
Anti-Sm . . . . .	24
Anti-SS-A 52 . . . . .	25
Anti-SS-A 60 . . . . .	26
Anti-SS-A . . . . .	27
Anti-SS-B . . . . .	28
Anti-ssDNA . . . . .	29
ENA-4-Profile . . . . .	30
ENA-6-Profile . . . . .	31
ENAcombi . . . . .	32
ENAscreen . . . . .	33
Rheumatoid Factor . . . . .	34

### Thrombosis Diagnostics

Anti-Annexin V . . . . .	35
Anti-beta-2-Glycoprotein I . . . . .	36
Anti-Cardiolipin . . . . .	37
Anti-Phosphatidic Acid . . . . .	38
Anti-Phosphatidyl Inositol . . . . .	39
Anti-Phosphatidyl Serine . . . . .	40
Anti-Phospholipid Screen . . . . .	41
Anti-Prothrombin . . . . .	42

### ANCA / Vasculitis Diagnostics

ANCAcombi . . . . .	43
ANCAscreen <i>hs (high sensitive)</i> . . . . .	44
Anti-BPI . . . . .	45
Anti-Cathepsin G . . . . .	46
Anti-Elastase . . . . .	47
Anti-GBM . . . . .	48
Anti-Lactoferrin . . . . .	49
Anti-Lysozyme . . . . .	50
Anti-MPO (pANCA) . . . . .	51
Anti-PR3 <i>hs</i> . . . . .	52

### Gastroenterology Diagnostics

AMA-M2 . . . . .	53
Anti-DGP . . . . .	54
Anti-Gliadin . . . . .	55
Anti-gp210 . . . . .	56
Anti-Intrinsic Factor . . . . .	57
Anti-LC1 . . . . .	58
Anti-LKM-1 . . . . .	59
Anti-Parietal Cell . . . . .	60
Anti-SLA . . . . .	61
Anti-Sp100 . . . . .	62
Anti-Tissue-Transglutaminase . . . . .	63
ASCA . . . . .	64

### Thyroid Diagnostics

Anti-TG . . . . .	65
Anti-TPO . . . . .	66
Thyroglobulin . . . . .	67

### Diabetes Diagnostics

Anti-Insulin . . . . .	68
------------------------	----

### Miscellaneous

beta-2-Microglobulin . . . . .	69
Ferritin . . . . .	70
Micro-Albumin . . . . .	71

## INFECTIOUS DISEASE DIAGNOSTICS

Anti-B. pertussis Toxin . . . . .	72
Anti-Borrelia . . . . .	73
Anti-Borrelia Liquor . . . . .	74
Anti-Chlamydia pneumoniae . . . . .	75
Anti-Chlamydia trachomatis . . . . .	76
Anti-EBV (EBNA-1) IgG . . . . .	77
Anti-EBV (VCA) . . . . .	78
Anti-EBV (VCA) IgG Liquor . . . . .	79
Anti-EBV (ZEBRA) IgM . . . . .	80
Anti-Helicobacter pylori . . . . .	81
Anti-Hepatitis E Virus . . . . .	82
Anti-HSV-1 . . . . .	83
Anti-HSV-1/2 . . . . .	84
Anti-HSV-1/2 IgG Liquor . . . . .	85
Anti-HSV-2 . . . . .	86
Anti-Measles Virus . . . . .	87
Anti-Measles Virus IgG Liquor . . . . .	88
Anti-Mumps Virus . . . . .	89
Anti-Mycoplasma pneumoniae . . . . .	90
Anti-Parvovirus B19 . . . . .	91
Anti-Rubella Virus IgG Liquor . . . . .	92
Anti-VZV . . . . .	93
Anti-VZV IgG Liquor . . . . .	94
Anti-Yersinia . . . . .	95

## STOOL DIAGNOSTICS

Alpha-1-Antitrypsin . . . . .	96
Calprotectin . . . . .	97
Lactoferrin . . . . .	98

## BONE AND MINERAL METABOLISM

25-OH Vitamin D <sub>3</sub> /D <sub>2</sub> . . . . .	99
25-OH Vitamin D <sub>3</sub> /D <sub>2</sub> Control Set . . . . .	100

## IMMUNOBLOTS

ANA-9-Line . . . . .	101
Myositis plus . . . . .	102
Nucleo-9-Line . . . . .	103
ANCA-3-Line . . . . .	104
Gastro-5-Line . . . . .	105
Liver-9-Line 2nd Generation . . . . .	106

## RAPID TEST

rheumachec® . . . . .	107
-----------------------	-----

## IMMUNOFLUORESCENCE

Anti-Endomysium Antibodies (AEA) . . . . .	108
Anti-nDNA Antibodies (nDNA) . . . . .	109
Anti-Neutrophil Cytoplasmic Antibodies (ANCA-Ethanol) . . . . .	110
Anti-Neutrophil Cytoplasmic Antibodies (ANCA-Formalin) . . . . .	111
Anti-Nuclear Antibodies HEp-2 (ANA-HEp-2) . . . . .	112
Anti-Smooth Muscle Antibodies (ASMA) . . . . .	113
Antoantibodies RL/RK/RS . . . . .	114

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— ANA Detect is an ELISA screening test for the qualitative determination of IgG autoantibodies against components of the cell nucleus and the cytoplasm in human serum or plasma.

— Some antinuclear antibodies (ANA) serve as markers for the diagnosis of rheumatic inflammatory connective tissue diseases like systemic lupus erythematosus (SLE), Sjogren's syndrome, Sharp syndrome (mixed connective tissue disease), polymyositis/dermatomyositis, and progressive systemic sclerosis (PSS), or CREST syndrome.

— Each individual well in the ANA Detect test is coated with a mixture of 26 human, recombinant, highly purified antigens. This combination guarantees high diagnostic sensitivity with excellent specificity and is the first choice for screening in cases of suspected systemic autoimmune disease.

— Patient samples that test positive with ANA Detect should subsequently be tested with the respective single ELISA tests or with ANAcombi from ORGENTEC.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 600	ORG 200-24
RANGE	qualitative	qualitative
CUT-OFF	Index > 1.2	Index > 1.2
CALIBRATION	CDC-Sera, WHO Wo/80	CDC-Sera, WHO Wo/80
COATING	SS-A 52 (Ro 52), SS-A 60 (Ro 60), SS-B (La), RNP/Sm, RNP-70, RNP-A, RNP-C, Sm-BB, Sm-D, Sm-E, Sm-F, Sm-G, Scl-70, Jo-1, dsDNA, ssDNA, polynucleosomes, mononucleosomes, histone proteins H1, H2A, H2B, H3, H4, Pm-Scl-100, centromere B	

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— ANAcombi is an ELISA test system for the qualitative determination of IgG antibodies against nine different components of the cell nucleus or the cytoplasm in human serum or plasma. It is used for the serological differentiation of inflammatory rheumatic connective tissue diseases.

— Antinuclear antibodies (ANA) are typically found in patients with systemic autoimmune diseases, such as systemic lupus erythematosus (SLE), Sjogren's syndrome, Sharp syndrome (mixed connective tissue disease), polymyositis/dermatomyositis, and progressive systemic scleroderma (PSS), or CREST syndrome. Some ANA are thus markers for these individual connective tissue diseases.

— Each row of twelve microplate wells is coated with one of the eight antigens. The isolated detection of individual autoantibodies with the ANAcombi test allows for the serological differentiation of these diseases.

SYSTEM	ELISA
TESTS	96 wells
CODE	ORG 539
RANGE	qualitative
CUT-OFF	Index > 1.2
CALIBRATION	CDC-Sera
COATING	SS-A 52 (Ro 52), SS-A 60 (Ro 60), SS-B (La), RNP-70, RNP/Sm, Sm, Jo-1, centromere B, Scl-70

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— ANAscreen is an ELISA test system for the qualitative measurement of IgG antibodies against nuclear antigens, SS-A 52 (Ro 52), SS-A 60 (Ro 60), SS-B (La), RNP-70, Sm, RNP/Sm, Scl-70, Centromere B, and Jo-1, in human serum or plasma samples.

— Antinuclear antibodies (ANA) are detected in patients with systemic autoimmune diseases such as systemic lupus erythematosus (SLE), Sjogren's syndrome, Sharp syndrome, polymyositis/dermatomyositis, and progressive systemic scleroderma (PSS). ANA thus partially act as makers for the individual connective tissue diseases.

— Each individual well of the ANAscreen test is coated with a mixture of nine human recombinant and native, highly purified antigens, allowing for sensitive autoantibody screening.

— After a positive result with this test, further differentiation of the autoantibodies should be carried out with the ANAcombi test or with the respective single antigen ELISA assays from ORGENTEC.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 538	ORG 238-24
RANGE	qualitative	qualitative
CUT-OFF	Index > 1.2	Index > 1.2
CALIBRATION	CDC-Sera, WHO Wo/80	CDC-Sera, WHO Wo/80
COATING	SS-A 52 (Ro 52), SS-A 60 (Ro 60), SS-B (La), RNP/Sm, RNP-70, Sm, Scl-70, Jo-1, centromere B	

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— Anti-alpha Fodrin IgA and Anti-alpha Fodrin IgG are ELISA-based, test system for the quantitative determination of IgA or IgG antibodies against alpha-fodrin in human serum or plasma.

— Alpha-fodrin is an intracellular, actin-binding, organ-specific protein of the cytoskeleton. During apoptosis alpha-fodrin is proteolytically cleaved. The resulting protein fragment accumulates in the salivary glands and may trigger a local autoimmune reaction.

— Antibodies against alpha-fodrin are specifically detected in serum samples from patients with primary or secondary Sjogren's syndrome. They can be detected very early, sometimes before appearance of SS-A (Ro) or SS-B (La) antibodies, which are important diagnostic markers for the disease.

— Sjogren's syndrome is an autoimmune disorder mainly affecting the salivary and lacrimal glands. Leading symptoms are dry eyes and mouth (sicca symptoms). Skin, mucous membranes, kidneys, blood vessels and other organs may also be affected.

SYSTEM	ELISA	Alegria	
TESTS	96 wells	24 Alegria® Test Strips	
CODE	ORG 642	ORG 242G-24	ORG 242A-24
RANGE	0–100 U/ml	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml	10 U/ml
CALIBRATION	internal	internal	internal
COATING	alpha-fodrin (human)		

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— Anti-C1q is an ELISA test system for the quantitative determination of IgG antibodies against the C1q complement factor in human serum or plasma.

— C1q is the first protein in the activation cascade of the classical complement pathway. Antibodies that bind to the collagen portion of C1q with high affinity have been detected in patients with a number of rheumatic disorders. In particular, anti-C1q antibodies have been found in patients with systemic lupus erythematosus (SLE) and patients with hypocomplementemic urticarial vasculitis syndrome (HUVS). For SLE patients there exists an unambiguous correlation between antibody titer, disease progression and kidney involvement. SLE patients with an increasing anti-C1q titer were observed to have an approximately 50 % higher risk of developing lupus nephritis. C1q antibodies often appear a few weeks, sometimes as much as half a year before onset of nephritis. Continuous monitoring of the antibody status with ORGENTEC's Anti-C1q test is thus recommended.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 549	ORG 249-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	C1q (human)	

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— Anti-CCP *hs (high sensitive)*® is an ELISA test system for the quantitative determination of IgG antibodies against cyclic citrullinated peptides (CCP) in human serum or plasma.

— Anti-CCP *hs (high sensitive)*® is based on selected and optimized peptide epitopes from the body's own MCV protein (modified citrullinated vimentin). Due to its tailored antigen profile, this ELISA attains the highest sensitivity while maintaining high specificity. It is thus positioned as a powerful serological test for the rapid and precise routine diagnosis of rheumatoid arthritis (RA), as well as for the diagnosis of early RA.

— In conjunction with rheumatoid factor, autoantibodies against citrullinated antigens (ACPA) have been included as determinate factors in the current ACR criteria for the classification of rheumatoid arthritis. Anti-CCP *hs (high sensitive)*® is an ACPA test that allows for the diagnosis of RA with comparable sensitivity and significantly higher specificity and prognostic power than conventional rheumatoid factor tests.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 601	ORG 301
RANGE	0–1000 U/ml	0–1000 U/ml
CUT-OFF	20 U/ml	20 U/ml
CALIBRATION	internal	internal
COATING	cyclic citrullinated vimentin peptides	

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— Anti-Centromere B is an ELISA test system for the quantitative determination of IgG antibodies against centromere B in human serum or plasma.

— The detection of autoantibodies against centromere protein B is of great diagnostic and prognostic importance in cases of progressive systemic scleroderma (PSS).

— Centromere protein B is a DNA-associated protein. Together with other proteins, it forms the kinetochores of the chromosomes. Centromere B is the target antigen in over 95 % of patients with scleroderma or Raynaud's phenomenon. Because autoantibodies against centromere proteins sometimes appear years before disease onset, detection of anti-centromere B can be a significant indicator for the development of scleroderma.

— CREST syndrome, a limited and localized form of scleroderma with a more favourable prognosis can be serologically detected at an early stage with use of the ORGENTEC Anti-Centromere B test.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 633	ORG 233-24
RANGE	0–300 U/ml	0–300 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	centromere B	

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— The anti-dsDNA ELISAs from ORGENTEC are quantitative assays for separate or combined determination of IgA, IgG, or IgM antibodies against double-stranded DNA in human serum or plasma samples.

— This assay is used for initial diagnosis when systemic lupus erythematosus is suspected and for differential diagnosis of other systemic autoimmune diseases. In addition to high titers of anti-nuclear antibodies (ANA), the detection of autoantibodies against dsDNA is an integral component of the ACR criteria for SLE diagnostics. Autoantibodies against DNA and histones bind to single-stranded (ssDNA) and double-stranded (dsDNA). Antibodies against dsDNA are predominantly detected during the active stages of the disease and serve as activity markers for SLE. They are a prognostic indicator of severe, multisystemic disease progression. The antibody concentration can also be used to monitor therapeutic success and to predict flares of the disease.

SYSTEM	ELISA				Alegria			
TESTS	96 wells				24 Alegria® Test Strips			
ORDER	ORG 604	ORG 604M	ORG 604A	ORG 604S	ORG 204G-24	ORG 204M-24	ORG 204A-24	ORG 204S-24
RANGE	0–200 IU/ml	0–200 U/ml	0–200 U/ml	0–200 U/ml	0–200 IU/ml	0–200 U/ml	0–200 U/ml	0–200 U/ml
CUT-OFF	20 IU/ml	20 U/ml	20 U/ml	20 U/ml	20 IU/ml	20 U/ml	20 U/ml	20 U/ml
CALIBRATION	WHO Wo/80	internal	internal	internal	WHO Wo/80	internal	internal	internal
COATING	double stranded DNA							

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— The Anti-Histone ELISA test allows for the quantitative determination of these antibodies of the IgG subtype in human serum or plasma. Up to 80 % of patients suffering from systemic lupus erythematosus (SLE) and nearly 95 % of patients with drug-induced lupus (DIL) have autoantibodies against histone proteins in their antibody profile.

— High titers of histone antibodies are almost without exception a sign of systemic lupus erythematosus or of drug-induced lupus, when other SLE markers are absent. When carrying out the Anti-Histone ELISA, parallel tests for dsDNA antibodies and Sm antibodies are recommended. They specifically indicate SLE.

— Aside from SLE and DIL, histone antibodies also occur in numerous other rheumatic diseases, including rheumatoid arthritis, juvenile chronic arthritis, and scleroderma..

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 507	ORG 207-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	40 U/ml	40 U/ml
CALIBRATION	WHO MRC 66/233	WHO MRC 66/233
COATING	histone proteins H1, H2A, H2B, H3, H4	

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— The Anti-Jo-1 ELISA is suited for the quantitative determination of IgG antibodies against Jo-1 in human serum or plasma. These antibodies are accepted, highly specific markers for idiopathic myositis.

— Histidyl-tRNA synthetase, which is one of the extractable nuclear antigens (ENA), is the antigen recognized by the Jo-1 autoantibodies.

— The diagnostic sensitivity of Anti-Jo-1 detection for idiopathic myositis in adults is as high as 30 %, for polymyositis it is up to 46 %. In cases of juvenile myositis, on the other hand, antibodies against Jo-1 occur less frequently. Their diagnostic specificity for autoimmune myositis is nearly 100 %. Anti-Jo-1 antibodies are absent in cases of neuromuscular disease without myositis.

— In myositis patients with a positive Jo-1 antibody test result, about 60 % of the cases involve antisynthetase syndrome with interstitial involvement of the lungs. Affected patients must expect a severe course of the disease with numerous flare-ups.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 513	ORG 213-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	histidyl-tRNA synthetase (Jo-1)	

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— Anti-MCV® is an ELISA test for the quantitative determination of IgG antibodies against mutated citrullinated vimentin (MCV) in human serum or plasma. It is a powerful serological test for rheumatoid arthritis, especially in early disease, even in patients who test negative for rheumatoid factor. Anti-MCV® antibodies may be detectable even years before the appearance of specific clinical symptoms.

— Autoantibodies against mutated citrullinated vimentin are highly specific and extremely sensitive markers for rheumatoid arthritis. Anti-MCV autoantibodies are an indicator of advancing damage to the joints and aggressive progression of the disease. In contrast to the detection of antibodies against synthetic cyclic citrullinated peptides, changes in the concentration of anti-MCV correspond to changes in clinical parameters. As well as allowing for an unambiguous diagnosis, the Anti-MCV® thus also makes it possible to monitor the response to treatment.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 548	ORG 248-24
RANGE	0–1000 U/ml	0–1000 U/ml
CUT-OFF	20 U/ml	20 U/ml
CALIBRATION	internal	internal
COATING	MCV (mutated citrullinated vimentin)	

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— Anti-Nucleosome is an ELISA test system for the quantitative determination of IgG antibodies against parts of the nucleosome, a structural component of chromatin, in human serum or plasma.

— This sensitive and specific serological test assists in the diagnosis of systemic lupus erythematosus (SLE) and drug-induced lupus (DIL). It is more sensitive than the test for dsDNA antibodies, and the anti-nucleosome antibodies can be detected earlier than anti-dsDNA antibodies. They thus provide a valuable indicator for the early diagnosis of emerging SLE in cases where SLE is suspected but dsDNA antibodies are not detected

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 528	ORG 228-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	20 U/ml	20 U/ml
CALIBRATION	internal	internal
COATING	nucleosomes (human)	

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— Anti-Rib-P is an ELISA test for the quantitative determination of IgG antibodies against ribosomal P proteins (Rib-P) in human serum or plasma. They recognize a common epitope of the three phosphoproteins P0, P1 and P2 of the 60S subunit of the ribosomal complex.

— Like anti-nuclear antibodies (ANA), anti-cytoplasmic antibodies play a role in rheumatic disease. Of these, the anti-ribosomal antibodies (ARA) are the best characterized. Autoantibodies against ribosomal P-proteins are highly specific markers for diagnostics of systemic lupus erythematosus, because they rarely occur in other autoimmune diseases.

— Anti-Rib-P detection is recommended in cases of suspected SLE when antibodies against Sm antigens and double-stranded DNA are absent. Use of the Anti-Rib-P ELISA from ORGENTEC is also indicated when lupus-induced psychosis is suspected.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 517	ORG 217-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	ribosomal P protein (Rib-P)	

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— Anti-RNP/Sm is an ELISA test system for the quantitative, determination of IgG antibodies against the RNP/Sm ribonucleoprotein complex in human serum or plasma.

— The RNP/Sm complex belongs to the extractable nuclear antigens (ENA) and is useful in the diagnosis of mixed connective tissue disorder (MCTD, Sharp syndrome) and related autoimmune diseases.

— The antigen is a multi-enzyme complex made of three U1-specific ribonucleoproteins (U1-A, U1-C, U1-70 kDa) and nine different Sm core proteins.

— Antibodies against the 70 kDa protein of this complex are a very specific marker for Sharp syndrome. Nearly all affected patients have autoantibodies against this target antigen; absence of RNP/Sm antibodies thus rules out this disease.

— The Sm proteins and the ribonucleoproteins A and C are recognized by antibodies that also occur in cases of mixed connective tissue disorder and systemic lupus erythematosus. However, SLE patients generally have lower titers of antibodies against the 70 kDa protein than MCTD patients.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 511	ORG 211-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	RNP/Sm complex	

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— Anti-RNP-70 is an ELISA test for the quantitative determination of IgG antibodies against RNP-70 in human serum or plasma.

— RNP-70, the 70 kDa protein of the U1-snRNP complex, is one of the extractable nuclear antigens (ENA). It is useful for the diagnosis of mixed connective tissue disease (MCTD, Sharp syndrome) and related autoimmune diseases.

— The diagnostic sensitivity of this marker for MCTD is 100 %; absence of RNP-70 autoantibodies rules out MCTD. Antibodies against RNP-70 are also found in up to 32 % of cases of systemic lupus erythematosus (SLE) and up to 10 % of systemic scleroderma cases.

— Patients whose only detectable antinuclear antibodies (ANA) are autoantibodies against the RNP-70 protein have the classical presentation of MCTD. If additional SLE-typical autoantibodies (anti-dsDNA, Anti-Sm, Anti-Histone) are found, transition into SLE is indicated.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 632	ORG 232-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	RNP-70 protein	

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— Anti-Scl-70 is an ELISA test for the quantitative determination of IgG antibodies against Scl-70 (DNA topoisomerase I) in human serum or plasma.

— With a specificity near 100 %, these antibodies are an accepted and highly specific marker for progressive systemic scleroderma; autoantibodies against topoisomerase I are extremely rare in both healthy individuals and those with non-autoimmune diseases.

— In addition to the serological diagnosis of systemic scleroderma, ORGENTEC's Anti-Scl-70 test is also useful in the differential diagnosis of Raynaud's phenomenon. Autoantibodies against topoisomerase I are also a prognostic marker for the further development of systemic sclerosis. Patients with Anti-Scl-70 antibodies have a higher tendency towards skin manifestations and involvement of the lungs, heart or kidney than patients who do not have these antibodies.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 512	ORG 212-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	Scl-70 protein (DNA topoisomerase I)	

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— Anti-Sm is an ELISA test for the quantitative determination of IgG antibodies against Sm proteins in human serum or plasma. The detection of autoantibodies against Sm proteins is one of the ACR criteria for the diagnosis of systemic lupus erythematosus (SLE).

— Because of the very high specificity (99 %) of this marker for SLE, the reliable detection of Sm autoantibodies suffices as a diagnostic criterion for the disease when fewer than three clinical ACR criteria have been met. However, a negative result does not rule out SLE.

— The detection of Sm antibodies serves as a prognostic marker for SLE. Obviously, there is a relationship between the appearance of Sm antibodies and severe organ manifestations of the disease.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 510	ORG 210-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	Sm proteins	

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— Anti-SS-A 52 is an ELISA test for the quantitative determination of IgG antibodies against SS-A 52 (Ro 52) in human serum or plasma. This test is used for the differential diagnosis and monitoring of systemic rheumatic inflammatory autoimmune diseases.

— Autoantibodies against SS-A 52 are primarily found in patients with collagenoses, especially those with Sjogren's syndrome and various forms of Lupus erythematosus. Autoantibodies against the SS-A 52 proteins are specific to Sjogren's syndrome, though they are less specific than SS-A 60 (Ro 60) antibodies. In cases of myositis, anti-SS-A 52 are found together with other myositis-specific antibodies such as Anti-SRP, Anti-PM-Scl, or antibodies against aminoacyl-tRNA synthetases, e.g. Jo-1. In these patients, there is no connection to Sjogren's syndrome.

— The involvement of the SS-A-52 autoantibodies in the pathogenesis of congenital heart block in lupus erythematosus seems likely.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 652	ORG 252-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	SS-A 52 (Ro 52)	

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— Anti-SS-A 60 is an ELISA test for the quantitative determination of IgG antibodies against SS-A 60 (Ro 60) in human serum or plasma. This test is useful for the differential diagnosis and monitoring of systemic rheumatic inflammatory autoimmune diseases.

— SS-A 60 is a component of the human cytoplasmic ribonucleoprotein complex. The proteins SS-A 60 and SS-A-52 bind to smaller RNAs and ribonucleoprotein particles in the cell nucleus and the cytoplasm. Autoantibodies against SS-A 60 are primarily found in cases of Sjögren's syndrome and lupus erythematosus. SS-A 60 antibodies frequently appear in conjunction with SS-B and SS-A 52 antibodies.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 660	ORG 260-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	SS-A 60 (Ro 60)	

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— Anti-SS-A is an ELISA test for the simultaneous quantitative determination of IgG antibodies against SS-A 52 (Ro 52) and SS-A 60 (Ro 60) in human serum or plasma. This test is useful for the differential diagnosis and monitoring of systemic rheumatic inflammatory autoimmune diseases.

— SS-A 52 and SS-A 60 are proteins that bind to smaller RNAs and ribonucleoprotein particles in the cell nucleus and the cytoplasm. They are thus a component of the human cytoplasmic ribonucleoprotein complex and are extractable nuclear antigens (ENA).

— Autoantibodies against SS-A 52 and SS-A 60 are predominantly found in cases of collagenoses, particularly Sjogren's syndrome, and various forms of lupus erythematosus.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 508	ORG 208-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	SS-A 52 (Ro 52), SS-A 60 (Ro 60)	

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— Anti-SS-B is an ELISA test for the quantitative determination of IgG antibodies against SS-B (La) in human serum or plasma.

— SS-B antigens are extractable nuclear antigens (ENA), and antibodies against SS-B are used for the differential diagnosis and monitoring of systemic inflammatory autoimmune diseases. Autoantibodies against the SS-B protein are found in cases of Sjogren's syndrome and various forms of lupus erythematosus. The diagnostic sensitivity of the marker for primary Sjogren's syndrome is 70 % and 50 % for secondary Sjogren's syndrome.

— SS-B antibodies are usually found together with anti-SS-A (Ro). In asymptomatic individuals who are suspected to have Sjogren's syndrome, SS-B antibodies can precede the full picture of the disease.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 509	ORG 209-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	SS-B (La)	

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— Anti-ssDNA is an ELISA test system for the quantitative, determination of IgG antibodies against single-stranded DNA in human serum or plasma. This assay is used for the serological detection of systemic lupus erythematosus (SLE) and drug-induced lupus (DIL).

— The antibody titers correlate with disease activity: In the serum of SLE patients, ssDNA antibodies are found in up to 87 % of patients during acute relapses of the disease but in only 43 % of patients during inactive phases.

— In conjunction with the detection of histone antibodies, the determination of ssDNA antibodies is helpful for the differential diagnosis of drug-induced lupus. ssDNA antibodies are found in over 50 % of these patients. Antibodies against ssDNA are additionally found in patients with collagenoses and juvenile chronic arthritis (JCA).

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 605	ORG 205-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	20 U/ml	20 U/ml
CALIBRATION	CDC sera	CDC sera
COATING	ssDNA (recombinant)	

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— ENA-4-Profile is an ELISA test system for the qualitative determination of IgG antibodies against four extractable nuclear antigens (ENA), SS-A (Ro 52, Ro 60), SS-B, Sm and RNP/Sm, in human serum or plasma. The test contributes to the diagnosis of inflammatory rheumatic connective tissue diseases.

— Antibodies against ENA are typically found in patients with lupus erythematosus (SLE), Sjogren's syndrome and mixed connective tissue disease. Some of these antibodies serve as markers for individual connective tissue diseases. Single detection of these antibodies with the ENA-4-Profile test thus makes it possible to differentiate these autoimmune disorders.

— The ENA-4-Profile test is based on a modular principle: Each kit contains three color-coded strips for each of the four antigens. Each strip has eight wells, each coated with an individual highly purified antigen.

SYSTEM	ELISA
TESTS	96 wells
CODE	ORG 544
RANGE	qualitative
CUT-OFF	Index > 1.2
CALIBRATION	CDC sera
COATING	RNP/Sm, Sm, SS-A (Ro 52, Ro 60), SS-B (La)

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— ENA-6-Profile is an ELISA test system for qualitative determination of IgG antibodies against six extractable nuclear antigens (ENA), SS-A (Ro 52, Ro 60), SS-B, Sm, RNP/Sm, Scl-70 and Jo-1 in human serum or plasma. The test aids to the diagnosis of systemic inflammatory rheumatic connective tissue diseases.

— Antibodies against ENA are typically found in patients with lupus erythematous (SLE), Sjogren's syndrome, Sharp syndrome, polymyositis/dermatomyositis, and progressive systemic scleroderma (PSS). Some of these antibodies serve as markers for individual connective tissue diseases. Single detection of these antibodies with the ENA-6-Profile test thus makes it possible to differentiate these diseases.

— The ENA-6-Profile test is based on a modular principle: Each kit contains two color-coded strips for each of the six antigens.

SYSTEM	ELISA
TESTS	96 wells
CODE	ORG 546
RANGE	qualitative
CUT-OFF	Index > 1.2
CALIBRATION	CDC sera
COATING	RNP/Sm, Sm, SS-A (Ro 52, Ro 60), SS-B (La), Scl-70, Jo-1

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— ENACombi is an ELISA test for the quantitative determination of IgG antibodies against SS-A 60, SS-A 52, SS-B, Sm, RNP/Sm, Scl-70, and Jo-1 in human serum or plasma. It is used in the diagnosis of inflammatory rheumatic connective tissue diseases.

— Autoantibodies against ENA can be detected in cases of systemic lupus erythematosus, Sjogren's syndrome, Sharp syndrome, polymyositis/dermatomyositis, and progressive systemic scleroderma. Some of these antibodies are disease markers and their isolated detection with ENACombi makes it possible to serologically differentiate between the individual connective tissue diseases.

— Each row of twelve microplate wells is coated with one of the six antigens (human recombinant or native affinity-purified). ENACombi thus allows for the determination of twelve differentiated patient profiles in one test.

SYSTEM	ELISA
TESTS	96 wells
CODE	ORG 514
RANGE	0–100 U/ml
CUT-OFF	25 U/ml
CALIBRATION	CDC sera
COATING	RNP/Sm, Sm, SS-A (Ro 52, Ro 60), SS-B (La), Scl-70, Jo-1

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— ENAscreen is an ELISA test system for the qualitative determination IgG antibodies against extractable nuclear antigens (ENA), SS-A (Ro), SS-B (La), Sm, RNP/Sm, Scl-70, Jo-1, in human serum or plasma. It is used in cases of suspected connective tissue disease.

— Besides DNA and histones, autoantibodies may also target other nuclear antigens. These nuclear antigens were named extractable nuclear antigens (ENA), as originally they were extracted from the nuclei with saline solution. Autoantibodies to Smith antigen (Sm) which is also considered to be highly specific for SLE were the first anti-ENA detected. Thereafter, further subtypes of ENA i.e. ribonucleoproteins (RNP), Scl-70, Jo-1, or the Sjogren's antigens SS-A and SS-B, Scl-70, and Jo-1 were identified.

— Each individual well of the ENAscreen test is coated with a mixture of the six highly purified antigens. For further differentiation, a positive result with this ELISA may then be followed up with ENAcombi or the ORGENTEC test kits for quantitative determination of the individual ENA specificities.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 506	ORG 206-24
RANGE	qualitative	qualitative
CUT-OFF	Index > 1.0	Index > 1.0
CALIBRATION	CDC sera	CDC sera
COATING	SS-A (Ro 52, Ro 60), SS-B (La), Sm, RNP/Sm, Scl-70, Jo-1	

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The ORGENTEC test kits Rheumatoid Factor IgG, IgM, IgA and Screen are ELISA test systems for the quantitative determination of rheumatoid factors (RF) of the individual immunoglobulin classes IgG, IgM, IgA, or any of the three subtypes simultaneously, in human serum or plasma.

Rheumatoid Factors are the best-known diagnostic markers for rheumatoid arthritis (RA), and presence of rheumatoid factors is an integral part of the current ACR criteria for classification of RA.

In established RA, high titer serum IgM RF indicate poor prognosis. The presence of either IgG or IgA RF in patients with long-standing RA may be a prognostic indicator of systemic manifestations.

For more extensive serological diagnosis, we recommend our Anti-MCV<sup>®</sup> ELISA or – especially for point-of-care diagnosis – the rheumachec<sup>®</sup> rapid test.

SYSTEM	ELISA				Alegria			
TESTS	96 wells				24 Alegria <sup>®</sup> Test Strips			
CODE	ORG 522G	ORG 522M	ORG 522A	ORG 522S	ORG 222G-24	ORG 222M-24	ORG 222A-24	ORG 222S-24
RANGE	0–500 U/ml	0–500 IU/ml	0–500 U/ml	0–500 U/ml	0–500 U/ml	0–500 IU/ml	0–500 U/ml	0–500 U/ml
CUT-OFF	20 U/ml	20 IU/ml	20 U/ml	25 U/ml	20 U/ml	20 IU/ml	20 U/ml	25 U/ml
CALIBRATION	internal	WHO 64/2	internal	internal	internal	WHO 64/2	internal	internal
COATING	Fc fragments (human)							

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— Anti-Annexin V IgG/IgM, Anti-Annexin V IgG, and Anti-Annexin V IgM are ELISA test systems for the quantitative measurement of IgG and IgM class autoantibodies against Annexin V in human serum or plasma. These assays contribute to the diagnosis of primary or secondary antiphospholipid syndrome (APS).

— Annexin V antibodies are found in patients with APS, SLE, or other collagenoses, as well as those with rheumatoid arthritis. It is also possible that these autoantibodies have pathogenic significance in pregnancy complications in APS patients.

— Anti-Annexin V antibodies may be detectable within the clinical framework of APS with a negative phospholipid antibody result. For comprehensive APS diagnostics, it is thus recommended to test for autoantibodies against Annexin V along with antibodies against cardiolipin, beta-2-glycoprotein I and prothrombin.

SYSTEM	ELISA	Alegria	
TESTS	96 wells	24 Alegria® Test Strips	
CODE	ORG 643	ORG 243G-24	ORG 243M-24
RANGE	0–100 U/ml	0–100 U/ml	0–100 U/ml
CUT-OFF	8 U/ml	8 U/ml	8 U/ml
CALIBRATION	internal	internal	internal
COATING	annexin V		

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— The ORGENTEC test kits Anti-beta-2-Glycoprotein I IgG, IgM, IgA and Screen are ELISA test systems for the quantitative determination of autoantibodies against beta-2-glycoprotein I of the individual immunoglobulin classes IgG, IgM, IgA, or any of the three subtypes simultaneously, in human serum or plasma. The assays contribute to the diagnosis of antiphospholipid syndrome (APS).

— The clinical symptoms of APS alone are not sufficiently specific to make a definitive diagnosis. Laboratory tests thus play an important role in the diagnosis of the disease. In patients with APS, autoantibodies are formed that bind to phospholipids like cardiolipin or to phospholipid-binding proteins like beta-2-glycoprotein. Detection of these autoantibodies is an integral part of the classification criteria issued by the International Society of Thrombosis and Hemostasis.

— beta-2-GPI antibodies are associated with the occurrence of arterial or venous thrombosis, as well as with recurring miscarriage. beta-2-GPI antibodies may be detected even when the cardiolipin test result is negative.

SYSTEM	ELISA			Alegria			
TESTS	96 wells			24 Alegria® Test Strips			
CODE	ORG 521	ORG 521A	ORG 521S	ORG 221G-24	ORG 221M-24	ORG 221A-24	ORG 221S-24
RANGE	0–100 U/ml	0–100 U/ml	0–90 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–90 U/ml
CUT-OFF	8 U/ml	8 U/ml	10 U/ml	8 U/ml	8 U/ml	8 U/ml	10 U/ml
CALIBRATION	Harris sera, IRP 97/656 (IgG), HCAL (IgG)/EY2C9 (IgM)			Harris sera, IRP 97/656 (IgG), HCAL (IgG)/EY2C9 (IgM)			
COATING	beta-2-glycoprotein I (beta-2-GPI)						

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— The ORGENTEC test kits Anti-Cardiolipin IgG, IgM, IgA and Screen are ELISA test systems for the quantitative determination of autoantibodies against cardiolipin of the immunoglobulin classes IgG, IgM, IgA, or of the three subtypes simultaneously, in human serum or plasma. The assays contribute to the diagnosis of antiphospholipid syndrome (APS).

— Regardless of their immunoglobulin class, cardiolipin autoantibodies are an important immunological classification criterion for systemic lupus erythematosus (SLE). In addition, the quantitative determination of anti-cardiolipin is very useful for monitoring SLE patients with secondary APS.

— ORGENTEC's anti-cardiolipin tests are all coated with cardiolipin and human beta-2-glycoprotein I as cofactor..

SYSTEM	ELISA			Alegria			
TESTS	96 wells			24 Alegria® Test Strips			
CODE	ORG 515	ORG 515A	ORG 515S	ORG 215G-24	ORG 215M-24	ORG 215A-24	ORG 215S-24
RANGE	0–120 GPL-U/ml 0–80 MPL-U/ml	0–120 APL-U/ml	0–90 U/ml	0–120 GPL-U/ml	0–80 MPL-U/ml	0–120 APL-U/ml	0–90 U/ml
CUT-OFF	10 GPL-U/ml 7 MPL-U/ml	10 APL-U/ml	10 U/ml	10 GPL-U/ml	7 MPL-U/ml	10 APL-U/ml	10 U/ml
CALIBRATION	Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)			Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)			
COATING	cardiolipin (cofactor: beta-2-GPI)						

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— Anti-Phosphatidic Acid IgG/IgM, Anti-Phosphatidic Acid IgG and Anti-Phosphatidic Acid IgM are ELISA test systems for the quantitative determination of IgG and IgM antibodies against phosphatidic acid in human serum or plasma.

— Phosphatidic acid antibodies are almost without exception detected together with cardiolipin and/or phosphatidylserine antibodies. Their presence is an indicator of antiphospholipid syndrome (APS) and is helpful in assessing the risk of thrombosis in patients with systemic lupus erythematosus (SLE) and lupus-like diseases.

— In addition to purified phosphatidic acid as antigen, the wells of this test are coated with human beta-2-glycoprotein as cofactor.

SYSTEM	ELISA	Alegria	
TESTS	96 wells	24 Alegria® Test Strips	
CODE	ORG 537	ORG 237G-24	ORG 237M-24
RANGE	0–100 GPL-U/ml, 0–100 MPL-U/ml	0–100 GPL-U/ml	0–100 MPL-U/ml
CUT-OFF	10 GPL-U/ml, 10 MPL-U/ml	10 GPL-U/ml	10 MPL-U/ml
CALIBRATION	Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)	Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)	
COATING	phosphatidic acid (cofactor: beta-2-GPI)		

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- Anti-Phosphatidyl Inositol IgG/IgM, Anti-Phosphatidyl Inositol IgG and Anti-Phosphatidyl Inositol IgM are ELISA test systems for the quantitative determination of IgG and IgM antibodies against phosphatidyl inositol in human serum or plasma.
- Autoantibodies against phosphatidyl inositol, a negatively charged phospholipid, generally appear together with cardiolipin antibodies and antibodies against phosphatidyl serine. However, anti-phosphatidyl inositol has also been detected by itself in patients with antiphospholipid syndrome (APS).
- In the absence of anti-cardiolipin and anti-phosphatidyl serine, the detection of phosphatidyl inositol autoantibodies is thus indicative of APS.
- In addition to the antigen phosphatidyl inositol, the wells of the ORGENTEC Anti-Phosphatidyl Inositol test are also coated with human beta-2-glycoprotein I as cofactor.

SYSTEM	ELISA	Alegria	
TESTS	96 wells	24 Alegria® Test Strips	
CODE	ORG 536	ORG 236G-24	ORG 236M-24
RANGE	0–100 GPL-U/ml, 0–100 MPL-U/ml	0–100 GPL-U/ml	0–100 MPL-U/ml
CUT-OFF	10 GPL-U/ml, 10 MPL-U/ml	10 GPL-U/ml	10 MPL-U/ml
CALIBRATION	Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)	Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)	
COATING	phosphatidyl inositol (cofactor: beta-2-GPI)		

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- Anti-Phosphatidyl Serine IgG/IgM, Anti-Phosphatidyl Serine IgG and Anti-Phosphatidyl Serine IgM are ELISA test systems for the quantitative determination of IgG and IgM antibodies against phosphatidyl serine in human serum or plasma.
- Phosphatidyl serine antibodies are associated with antiphospholipid syndrome.
- Anti-cardiolipin and antibodies against phosphatidyl serine do not always appear together. When APS is suspected and anti-cardiolipin is absent, a further test for autoantibodies against phosphatidyl serine is particularly recommended.
- In addition to the antigen phosphatidyl serine, the wells of the ORGENTEC Anti-Phosphatidyl Serine tests are coated with human beta-2-glycoprotein I as cofactor.

SYSTEM	ELISA	Alegria	
TESTS	96 wells	24 Alegria® Test Strips	
CODE	ORG 535	ORG 235G-24	ORG 235M-24
RANGE	0–100 GPL-U/ml, 0–100 MPL-U/ml	0–100 GPL-U/ml	0–100 MPL-U/ml
CUT-OFF	10 GPL-U/ml, 10 MPL-U/ml	10 GPL-U/ml	10 MPL-U/ml
CALIBRATION	Harris sera, IRP 97/656 (IgG), HCAL (IgG)/EY2C9 (IgM)	Harris sera, IRP 97/656 (IgG), HCAL (IgG)/EY2C9 (IgM)	
COATING	phosphatidyl serine (cofactor: beta-2-GPI)		

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— Anti-Phospholipid Screen IgG/IgM, IgG and IgM are ELISA screening tests for the quantitative determination of IgG and IgM antibodies against cardiolipin, phosphatidyl serine, phosphatidyl inositol and phosphatidic acid in human serum or plasma. The wells of the microplate are coated with a mixture of the four phospholipids and beta-2-glycoprotein I as cofactor.

— This test is an outstanding tool for the serological diagnosis of primary and secondary antiphospholipid syndrome. The combined detection of four “phospholipid antibodies” lends this test particularly high clinical sensitivity.

— ORGENTEC’s test systems for the isolated detection of the individual phospholipid autoantibodies, and the Anti-beta-2-glycoprotein-I test can subsequently be used for the differentiated quantitative determination of these autoantibody species.

SYSTEM	ELISA	Alegria	
TESTS	96 wells	24 Alegria® Test Strips	
CODE	ORG 529	ORG 229G-24	ORG 229M-24
RANGE	0–100 GPL-U/ml, 0–100 MPL-U/ml	0–100 GPL-U/ml	0–100 MPL-U/ml
CUT-OFF	10 GPL-U/ml, 10 MPL-U/ml	10 GPL-U/ml	10 MPL-U/ml
CALIBRATION	Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)	Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)	
COATING	cardiolipin, phosphatidic acid, phosphatidyl inositol, phosphatidyl serine (cofactor: beta-2-GPI)		

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The ORGENTEC test kits Anti-Prothrombin IgG, IgM, IgA and Screen are ELISA test systems for the quantitative determination of autoantibodies against prothrombin of the immunoglobulin classes IgG, IgM, IgA, or of the three subtypes simultaneously, in human serum or plasma. The assays contribute to the diagnosis of antiphospholipid syndrome (APS).

Antiphospholipid syndrome (APS, Hughes syndrome) is a systemic autoimmune disease that may cause thromboses, recurring miscarriage, and intrauterine fetal death. In addition to these clinical symptoms, specific autoantibodies can be detected in blood samples from APS patients. These antibodies bind to phospholipids like cardiolipin or phospholipid-binding proteins like beta-2-glycoprotein.

The detection of prothrombin autoantibodies supports the diagnosis of antiphospholipid syndrome (APS) especially when cardiolipin antibodies are not detected. In high-risk groups for secondary APS such as SLE and scleroderma patients, these tests make it possible to assess the risk of thrombosis.

SYSTEM	ELISA			Alegria			
TESTS	96 wells			24 Alegria® Test Strips			
CODE	ORG 541	ORG 541A	ORG 541S	ORG 241G-24	ORG 241M-24	ORG 241A-24	ORG 241S-24
RANGE	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml	20 U/ml	10 U/ml	10 U/ml	10 U/ml	20 U/ml
CALIBRATION	Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)			Harris sera, IRP 97/656 (IgG), HCAL (IgG) / EY2C9 (IgM)			
COATING	prothrombin						

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— ANCAcombi is an ELISA test system for the qualitative determination of anti-neutrophil cytoplasmic antibodies (ANCA) in human serum or plasma. It is used for the serological detection of autoimmune systemic vasculitis.

— Antineutrophil cytoplasmic antibodies (ANCA) are directed against cytoplasmic antigens from human neutrophil granulocytes. They are primarily detected in patients with inflammatory disease of the vascular system. The ANCAcombi assay detects ANCA against seven antigens in a single test: proteinase 3 (PR3), myeloperoxidase (MPO), bactericidal permeability-increasing protein (BPI), elastase, cathepsin G, lysozyme, and lactoferrin.

— Each row of twelve microplate wells is coated with one of the seven antigens. An additional row is used for test calibration. Each ANCAcombi test kit thus allows to test twelve serum samples for the individual ANCA antigens.

— The corresponding ORGENTEC single ELISA tests can subsequently be used to quantitatively identify the individual ANCA antibodies.

<b>SYSTEM</b>	<b>ELISA</b>
<b>TESTS</b>	<b>96 wells</b>
<b>CODE</b>	<b>ORG 530</b>
RANGE	qualitative
CUT-OFF	Index $\geq$ 1.0
CALIBRATION	internal
COATING	proteinase 3 (PR3), myeloperoxidase (MPO), bactericidal permeability-increasing protein (BPI), cathepsin G, elastase, lactoferrin, lysozyme

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— ANCAscreen *hs* (high sensitive) is an ELISA test system for the qualitative determination of IgG antibodies against proteinase 3 (PR3) and myeloperoxidase (MPO) in human serum or plasma.

— Anti-PR3 and Anti-MPO are well defined serologic markers for ANCA-associated vasculitides (AAV).

— PR3 is the most frequent component of cANCA and the landmark autoantigen in granulomatosis with polyangiitis (GPA, formerly named Wegener's granulomatosis) with a clinical specificity of more than 95 %. MPO, the main target antigen of pANCA, is present in 70 % of patients with microscopic polyangiitis (MPA) and differentiates MPA from other autoimmune diseases.

— The wells of the microplate of the ANCAscreen *hs* (high sensitive) test are coated with a mixture of highly purified human MPO and PR3. For coating with PR3 we use the same innovative anchor technology, which we have developed and implemented in the highly sensitive single antigen ELISA Anti-PR3 *hs*. This high diagnostic performance has thus been integrated in the screening test.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 689	ORG 389
RANGE	qualitative	qualitative
CUT-OFF	Index ≥ 1.0	Index ≥ 1.0
CALIBRATION	internal	internal
COATING	proteinase 3 (PR3), myeloperoxidase (MPO)	

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— Anti-BPI is an ELISA-based test system for the quantitative determination of IgG antibodies against bactericidal permeability-increasing protein (BPI) in human serum or plasma.

— Antineutrophil cytoplasmic antibodies (cANCA) against bactericidal permeability-increasing protein (BPI), in addition to PR3 or MPO antibodies, are important serological markers in the diagnosis of a number of inflammatory, rheumatic and gastrointestinal diseases.

— Antibodies against BPI are primarily found in patients with cystic fibrosis and in patients with chronic inflammatory intestinal conditions like Crohn's disease and ulcerative colitis. Anti-PR3-negative serum samples that display a typical cANCA pattern in their immunofluorescence test should be tested for BPI antibodies.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 523	ORG 223-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	bactericidal permeability-increasing protein (BPI)	

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— Anti-Cathepsin G is an ELISA-based test system for the quantitative determination of IgG antibodies against cathepsin G in human serum or plasma. It is used for the differentiation of a positive pANCA immunofluorescence pattern. Cathepsin G antibodies are found in cases of systemic vasculitis, inflammatory rheumatic diseases (SLE, Sjogren's syndrome, Felty's syndrome), autoimmune hepatitis, and chronic inflammatory bowel diseases.

— Cathepsins are intracellular proteases of the lysosome. Their hydrolytic activity significantly contributes to the destruction of the bone matrix. Autoantibodies against cathepsin G give a perinuclear (pANCA) or atypical (xANCA) fluorescence pattern in indirect fluorescence on ethanol-fixed granulocytes.

— Patient sera that display such an immunofluorescence pattern but give a negative Anti-MPO ELISA result should be tested for cathepsin G antibodies, among others.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 525	ORG 225-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	cathepsin G	

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— Anti-Elastase is an ELISA test for the quantitative determination of IgG antibodies against elastase in human serum or plasma. These autoantibodies are particularly common in patients with chronic inflammatory bowel diseases and autoimmune hepatitis.

— Elastase is a serine protease. Its proteolytic effect plays a role in the destruction of tissue that is associated with rheumatoid arthritis and pulmonary emphysema.

— Like autoantibodies against myeloperoxidase (MPO), elastase autoantibodies are antineutrophil cytoplasmic antibodies, which evince a perinuclear fluorescence pattern (pANCA) in indirect immunofluorescence. Serum samples that display such a pANCA pattern but give a negative Anti-MPO ELISA result should be tested for autoantibodies against elastase, as well as antibodies against cathepsin G, lysozyme, and lactoferrin.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 524	ORG 224-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	elastase	

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— Anti-GBM is an ELISA test for the quantitative determination of IgG antibodies against components of the glomerular basement membrane (GBM) in human serum or plasma.

— These antibodies bind to the C-terminal end of the alpha-3 chain of type IV collagen (NC1 alpha-3 IV, the so called "Goodpasture antigen").

— The detection of these autoantibodies is unequivocal serological evidence for the diagnosis of Goodpasture syndrome, an autoimmune disease of the kidneys. The wells of the Anti-GBM test are coated with the highly purified antigen of the glomerular basement membrane. This guarantees high sensitivity and specificity of the ORGENTEC ELISA.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 550	ORG 250-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	20 U/ml	20 U/ml
CALIBRATION	internal	internal
COATING	glomerular basement membrane (GBM)	

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- Anti-Lactoferrin is an ELISA test system for the quantitative determination of IgG antibodies against lactoferrin in human serum or plasma.
- Autoantibodies against lactoferrin are classed as pANCA. They are often found in cases of inflammatory rheumatic and gastrointestinal autoimmune diseases such as rheumatoid vasculitis, or ulcerative colitis. High titers of lactoferrin are also found in autoimmune liver diseases.
- The Anti-Lactoferrin ELISA from ORGENTEC facilitates the quantitative serological detection of these autoantibodies, and contributes to the clarification of negative proteinase-3 and myeloperoxidase ANCA results.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 527	ORG 227-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	lactoferrin	

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— Anti-Lysozyme is an ELISA test for the quantitative determination of IgG antibodies against lysozyme in human serum or plasma.

— This test for antineutrophil cytoplasmic antibodies (ANCA) against lysozyme allows cases involving a positive immunofluorescence test for ANCA, but negative direct tests for proteinase 3 (PR3) and myeloperoxidase (MPO) antibodies to be cleared up quickly and reliably.

— Lysozyme is a glucosidase localized in the specific granules of neutrophils. In indirect immunofluorescence tests, lysozyme antibodies give a typical pANCA pattern. Their detection indicates rheumatoid vasculitis or Sjogren's syndrome; they may also occur in cases of SLE and ulcerative colitis.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 226-24</b>
RANGE	0–100 U/ml
CUT-OFF	10 U/ml
CALIBRATION	internal
COATING	lysozyme

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- Anti-MPO is an ELISA test for the quantitative determination of IgG antibodies against myeloperoxidase (MPO) in human serum or plasma. The anti-MPO assay is an important element in the diagnosis of microscopic polyangiitis (MPA).
- MPO antibodies are also found in patients with eosinophilic granulomatosis with polyangiitis (EGPA, Churg-Strauss syndrome) and polyarteritis nodosa. They can be used to serologically differentiate these diseases from granulomatosis with polyangiitis (GPA, formerly Wegener's disease).
- Antineutrophil cytoplasmic antibodies (ANCA) are a group of autoantibodies directed against cytoplasmic components of the neutrophils and monocytes. Current recommendations for ANCA detection include initial testing by indirect immunofluorescence (IF), followed by an antigen-specific ELISA such as Anti-MPO.
- MPO is a lysosomal enzyme of neutrophil granulocytes and monocytes. Antibodies against MPO exclusively recognize native conformational epitopes. The wells of the Anti-MPO ELISA are coated with highly purified human myeloperoxidase.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 519	ORG 219-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	5 U/ml	5 U/ml
CALIBRATION	internal	internal
COATING	myeloperoxidase	

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— Anti-PR3 *hs* (*high sensitive*) (ORG 618, ORG 318-24), Anti-PR3 (ORG 218-24) and Anti-PR3 (cANCA) (ORG 518) are ELISA tests for the quantitative determination of IgG antibodies against proteinase 3 (PR3) in human serum or plasma.

— The Anti-PR3 assay is an important element in the diagnosis of granulomatosis with polyangiitis (formerly Wegener's disease).

— Based on a special coating technology, Anti-PR3 *hs* now represents an innovative ELISA whose extraordinarily high sensitivity and specificity provide outstanding diagnostic significance in ANCA serology. In addition, the results of Anti-PR3 *hs* correlate well with those of indirect immunofluorescence (IIF).

— Current guidelines recommend an antigen-specific ELISA for the detection of autoantibodies against PR3, in combination with and IIF test if needed.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 618	ORG 318-24
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	proteinase 3	

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— AMA-M2 is an ELISA test for the quantitative determination of IgG antibodies against the mitochondrial antigen of the subtype M2 in human serum or plasma.

— Anti-mitochondrial antibodies (AMA) are a heterogeneous group of autoantibodies directed against various proteins of the outer and inner mitochondrial membrane. Autoantibodies of the M2 subtype bind to the E2 subunits of the pyruvate dehydrogenase complex (PDC-E2), the branched chain-2-oxoacid dehydrogenase complex (BCOADC-E2), and the 2-oxoglutarate dehydrogenase complex E (OGDC-E2).

— AMA-M2 are highly sensitive and specific. They represent a serological hallmark of primary biliary cirrhosis. In addition to early detection and differential diagnosis, determination of the AMA profile allows for an immunological and prognostic classification of primary biliary cirrhosis.

— In patients with other autoimmune diseases, especially collagenoses, detection of AMA-M2 makes it possible to spot a developing or associated primary biliary cirrhosis.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 516	ORG 216-24
RANGE	0–200 IU/ml	0–200 IU/ml
CUT-OFF	10 IU/ml	10 IU/ml
CALIBRATION	WHO 67/183	WHO 67/183
COATING	M2 proteins (PDC-E2, BCOADC-E2, OGDC-E2)	

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Anti-DGP IgA is an ELISA test for the quantitative determination of IgA and/or IgG antibodies against deamidated gliadin peptide epitopes (DGP) in human serum or plasma. The assays Anti-DGP IgA, Anti-DGP IgG, and Anti-DGP Screen are used to support a diagnosis of celiac disease.

After resorption, gliadin from food is deamidated by tissue transglutaminase (tTG) in the gastric mucosa. Deamidated gliadin peptides bind to the HLA molecules of antigen-presenting cells. The complex is recognized by corresponding T-helper cells. This results in an extensive immune response that includes formation of antibodies against tTG and DGP.

Because a substantial percentage of celiac patients are deficient in IgA antibodies, sensitive screening of individuals at risk always also includes testing for IgG antibodies. Anti-DGP Screen detects antibodies of both subclasses, IgA and IgG, in a single test step, and thus also identifies celiac disease patients with IgA deficiency

SYSTEM	ELISA			Alegria		
TESTS	96 wells			24 Alegria® Test Strips		
CODE	ORG 551G	ORG 551A	ORG 551S	ORG 251G	ORG 251A	ORG 251S
RANGE	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml	10 U/ml	10 U/ml	10 U/ml	10 U/ml
CALIBRATION	internal	internal	internal	internal	internal	internal
COATING	deamidated gliadin proteins					

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— The Anti-Gliadin ELISA tests from ORGENTEC are used for the quantitative determination of IgA and/or IgG antibodies against gliadin in human serum or plasma. The assays Anti-Gliadin IgA, Anti-Gliadin IgG, and Anti-Gliadin Screen support a diagnosis of celiac disease. They are helpful in monitoring maintenance of a gluten-free diet, which is the only effective therapy.

— The protein gliadin is a component of gluten, a mixture of proteins found in many types of grain. A permanent intolerance to gluten is the cause of celiac disease, which causes villous atrophy and an associated reduction in the absorption of nutrients.

— Anti-Gliadin titers correlate with the inflammation and with the condition of the mucous membrane of the small intestine. A gluten-free diet results in decrease in the antibody titer and thus leads to a remission of the disease.

SYSTEM	ELISA			Alegria		
TESTS	96 wells			24 Alegria® Test Strips		
CODE	ORG 534G	ORG 534A	ORG 534S	ORG 234G-24	ORG 234A-24	ORG 234S-24
RANGE	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml	0–100 U/ml
CUT-OFF	12 U/ml	12 U/ml	15 U/ml	12 U/ml	12 U/ml	15 U/ml
CALIBRATION	internal	internal	internal	internal	internal	internal
COATING	gliadin					

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— Anti-gp210 is an ELISA-based, automated test system for the quantitative determination of IgG antibodies against nuclear glycoprotein 210 (gp210) in human serum or plasma. It is used in suspected cases of primary biliary cirrhosis (PBC).

— Glycoprotein 210 is a component of the nuclear pore complex. Autoantibodies against gp210 are highly specific to primary biliary cirrhosis and indicate an unfavorable prognosis for the PBC patient.

— AMA-M2 detection is the primary diagnostic marker for PBC; however, AMA are absent in some patients. In such cases, complementary detection of antibodies against gp210, Sp100, or centromere can be diagnostically helpful.

— Antibodies against gp210 are very rarely also found in patients with autoimmune hepatitis, rheumatoid arthritis, polymyositis, or Sjogren's syndrome.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 256</b>
RANGE	0–100 U/ml
CUT-OFF	10 U/ml
CALIBRATION	internal
COATING	nuclear pore antigen gp210

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— Anti-Intrinsic Factor is an ELISA test for the quantitative determination of IgG antibodies against intrinsic factor in human serum or plasma.

— The detection of these antibodies aids in the diagnosis of pernicious anemia. They are found in 50–70 % of affected patients and are highly specific to this disease.

— Intrinsic factor is a protein formed by parietal cells in the stomach mucosa. The factor plays a crucial role in the absorption of vitamin B12 in the small intestine. Pernicious anemia causes damage to the gastric mucosa. This inhibits the synthesis of intrinsic factor and results in a vitamin B12 deficiency. These patients also generate autoantibodies against parietal cells and intrinsic factor.

— The Anti-Intrinsic Factor test from ORGENTEC uses the human recombinant intrinsic factor as antigen and thus guarantees maximum efficiency for the differential diagnosis of vitamin B12 deficiency.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 647	ORG 247-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	6 U/ml	6 U/ml
CALIBRATION	internal	internal
COATING	intrinsic factor (human recombinant)	

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— Anti-LC1 is an ELISA-based, automated, in-vitro test system for the quantitative determination of IgG antibodies against liver cytosolic protein type 1 (LC1) in human serum or plasma. It is used to support the differential diagnosis of liver diseases, particularly in cases of suspected type 2 autoimmune hepatitis.

— The detection of autoantibodies plays a central role in the diagnosis of inflammatory liver diseases namely autoimmune hepatitis (AIH), primary biliary cirrhosis and sclerosing cholangitis in adults and children. Anti-smooth muscle antibodies (ASMA) and anti-nuclear antibodies (ANA) define AIH type 1. Liver-kidney microsome antibodies (anti-LKM-1) and liver cytosolic protein type 1 antibodies (anti-LC1) may be present in AIH type 2.

— Anti-LC1 may be the single antibody in up to one third of the patients with AIH-2. Anti-LC1 antibodies have also been found in 14 % anti-LKM-1 antibody positive patients suffering from chronic hepatitis C virus infection.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 258</b>
RANGE	0–100 U/ml
CUT-OFF	10 U/ml
CALIBRATION	internal
COATING	LC1 protein (recombinant)

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— Anti-LKM-1 is an ELISA-based test system for the quantitative determination of IgG antibodies against liver kidney microsome type 1 (LKM-1) in human serum or plasma. It is used to support the differential diagnosis of liver diseases, particularly in cases of suspected type 2 autoimmune hepatitis.

— Anti-smooth muscle antibodies (ASMA) and anti-nuclear antibodies (ANA) are specific markers for AIH 1 anti-LKM-1 and antibodies to liver cytosol type 1 (anti-LC1) indicate type 2 AIH.

— ANA, ASMA, anti-LKM-1 and AMA-M2 should be determined in all patients suggestive of autoimmune liver disease. Monitoring of autoantibodies may be useful as disappearance or sharp decrease of ANA, ASMA and anti-LKM-1 can be an indicator of response to immunosuppressive treatment.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 253</b>
RANGE	0–100 U/ml
CUT-OFF	10 U/ml
CALIBRATION	internal
COATING	LKM-1 protein (recombinant)

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— Anti-Parietal Cell is an ELISA test for the quantitative determination of IgG antibodies against the H<sup>+</sup>/K<sup>+</sup>-ATPase of parietal cells in human serum or plasma.

— This enzyme transports hydrogen ions into the stomach lumen, regulating the acidic pH value of the gastric juice.

— ELISA test systems for parietal cell antigen are distinguished by their high sensitivity and specificity. Circulating autoantibodies against H<sup>+</sup>/K<sup>+</sup>-ATPase are found in 80–90 % of patients with pernicious anemia; however, they are also detected in 2–5 % of healthy individuals. In general, parietal cell autoantibodies are found more frequently in older people than in younger ones.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 531	ORG 231-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	H <sup>+</sup> /K <sup>+</sup> -ATPase (human recombinant)	

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
— Anti-SLA is an ELISA-based test system for the quantitative determination of IgG antibodies against soluble liver antigen (SLA) in human serum or plasma. It is useful in the diagnosis of liver disease, particularly when autoimmune hepatitis type 1 (AIH-1) is suspected.

— Antibodies against SLA (soluble liver antigen; synonym: liver pancreas antigen, LP) are the most specific markers for AIH-1. Anti-SLA is often found in conjunction with ASMA or ANA. If these typical AIH antibodies are absent despite clinical presentation of AIH, the isolated detection of SLA antibodies takes on special diagnostic significance. In addition, SLA antibodies are indicative of a more severe progression of autoimmune hepatitis.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 257</b>
RANGE	0–100 U/ml
CUT-OFF	10 U/ml
CALIBRATION	internal
COATING	SLA/LP protein

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— Anti-Sp100 is an ELISA-based test system for the quantitative determination of IgG antibodies against soluble nuclear protein Sp100 in human serum or plasma. It is used when primary biliary cirrhosis is suspected.

— Primary biliary cirrhosis (PBC) is primarily indicated by the detection of anti-mitochondrial antibodies of the M2 subtype (AMA-M2 antibodies). However, in some PBC patients, these AMA are absent. In such cases, SP100 antibody detection is a very valuable tool in the diagnostic kit.

— At 95 %, the diagnostic specificity of Sp100 antibodies toward PBC is very high; they are very rarely detected in patients with autoimmune hepatitis (AIH) or primary sclerosing cholangitis (PSC). The simultaneous detection of AMA-M2 and Sp100 antibodies is highly specific to PBC.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 254</b>
RANGE	0–100 U/ml
CUT-OFF	10 U/ml
CALIBRATION	internal
COATING	Sp100 protein (recombinant)

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The Anti-Tissue-Transglutaminase ELISA tests from ORGENTEC are used for quantitative determination of IgA and/or IgG antibodies against Tissue-Transglutaminase (tTG) in human serum or plasma. The assays Anti-Tissue-Transglutaminase IgA, Anti-Tissue-Transglutaminase IgG, and Anti-Tissue-Transglutaminase Screen support a diagnosis of celiac disease.

tTG belongs to a family of calcium-dependent enzymes that catalyze the transfer of protein-bound glutamine residues to primary amines. The enzyme is released by inflamed cells and is the primary antigen for autoantibodies of celiac disease patients.

Due to its high reliability, sensitivity, and specificity, the Anti-Tissue-Transglutaminase assay is an ideal non-invasive test for screening patients with clinical symptoms. International Guidelines, e.g. ESPGAN, NICE recommend the test for serological diagnosis of celiac disease.

SYSTEM	ELISA			Alegria		
TESTS	96 wells			24 Alegria® Test Strips		
CODE	ORG 540G	ORG 540A	ORG 540S	ORG 240G-24	ORG 240A-24	ORG 240S-24
RANGE	0–200 U/ml	0–200 U/ml	0–200 U/ml	0–200 U/ml	0–200 U/ml	0–200 U/ml
CUT-OFF	10 U/ml	10 U/ml	15 U/ml	10 U/ml	10 U/ml	15 U/ml
CALIBRATION	internal	internal	internal	internal	internal	internal
COATING	tissue transglutaminase (tTG) (human recombinant)					

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— ASCA IgG/IgA, ASCA IgG and ASCA IgA are ELISA based tests for the quantitative determination of IgG and IgA antibodies against *Saccharomyces cerevisiae* (ASCA) in human serum or plasma. These antibodies are directed against mannan in the cell membrane of *S. cerevisiae*.

— ASCA detection facilitates the differential diagnosis of the inflammatory bowel diseases, Crohn's disease and ulcerative colitis, particularly in the early stages of disease.

— The serum ASCA concentration is significantly higher in patients with Crohn's disease than in those with ulcerative colitis. Studies show a sensitivity of 75 % toward IgG-ASCA and 60 % toward IgA-ASCA in patients with Crohn's disease. In contrast IgG titers are elevated in only 5 % of patients with ulcerative colitis, and IgA titers in only 7 %.

— An additional test for perinuclear anti-cytoplasmic antibodies (pANCA) improves differential diagnosis: the prevalence of pANCA in ulcerative colitis is 50–90 %, while in Crohn's disease it is only 10–20 %.

SYSTEM	ELISA	Alegria	
TESTS	96 wells	24 Alegria® Test Strips	
CODE	ORG 545	ORG 245A-24	ORG 245G-24
RANGE	0–100 U/ml	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml	10 U/ml
CALIBRATION	internal	internal	internal
COATING	mannan ( <i>Saccharomyces cerevisiae</i> )		

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- Anti-TG is an ELISA test for the quantitative determination of IgG antibodies against thyroglobulin in human serum or plasma.
- The protein thyroglobulin is responsible for the production and storage of thyroid hormones and is a strong autoantigen. B-lymphocytes with membrane-bound IgM antibodies against thyroglobulin are even detected in the blood of healthy individuals.
- High titers of IgG antibodies against thyroglobulin found in conjunction with antibodies against thyroid peroxidase (Anti-TPO) are markers for Hashimoto-type chronic autoimmune thyroiditis (Hashimoto's disease). Thyroglobulin antibodies are found in 50–60 % of patients with this condition.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 502	ORG 202-24
RANGE	0–9000 IU/ml	0–9000 IU/ml
CUT-OFF	150 IU/ml	150 IU/ml
CALIBRATION	WHO 65/93	WHO 65/93
COATING	thyroglobulin (human)	

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— Anti-TPO is an ELISA test for the quantitative determination of IgG antibodies against thyroid peroxidase (TPO) in human serum or plasma.

— Thyroid peroxidase is important for the production of thyroid hormones in that it regulates the oxidative iodination of the tyrosine residues in thyroglobulin. If TPO is inhibited by binding of specific antibodies, less thyroid hormones are produced, which leads to hypothyroidism. Because of the close functional relationship between the two enzymes, antibodies against TPO and thyroglobulin are often found together.

— Elevated TPO antibody titers are found in 90 % of patients with autoimmune thyroid disease (Hashimoto's thyroiditis). TPO antibodies can also be detected in patients with autoimmune hyperthyroidism (Grave's disease). In some cases, trace concentrations of TPO antibodies are found in healthy individuals, particularly those of advanced age.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 503	ORG 203-24
RANGE	0–3000 IU/ml	0–3000 IU/ml
CUT-OFF	75 IU/ml	75 IU/ml
CALIBRATION	MRC 66/387	MRC 66/387
COATING	thyroid peroxidase (human)	

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- Thyroglobulin is an ELISA test for the quantitative determination of thyroglobulin (hTG) in human serum or plasma. It is used for the diagnosis and monitoring of pathological changes in the thyroid gland.
- Thyroglobulin is produced in the thyroid gland and is responsible for the synthesis of thyroid hormones T3 and T4. Small amounts of thyroglobulin also continually enter the blood stream.
- Elevated thyroglobulin levels are found in patients with pathological changes to the thyroid gland, for example hyperthyroidism, euthyroid goiter, thyroiditis, or thyroid carcinoma.
- The thyroglobulin ELISA is primarily used for postoperative monitoring of patients with differentiated thyroid carcinoma. After a successful thyroidectomy and targeted radiation therapy to remove residual thyroid tissue, the hTG concentration should soon sink below the detection limit after an initial temporary increase.

<b>SYSTEM</b>	<b>ELISA</b>
<b>TESTS</b>	<b>96 wells</b>
<b>CODE</b>	<b>ORG 5TG</b>
RANGE	0–300 ng/ml
CALIBRATION	internal
COATING	anti-thyroglobulin antibodies

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— Anti-Insulin is an ELISA test for the quantitative determination of IgG antibodies against human, bovine and pig insulin in human serum or plasma.

— Type I diabetes mellitus is an autoimmune disease triggered by autoantibodies against pancreatic islet cells and against insulin itself. Autoantibodies against insulin can appear in the prediabetic phase of type I diabetes mellitus or once the condition manifests. These autoantibodies are found in 50 to 70 % of children with diabetes mellitus and 20 to 30 % of affected adults.

— In rare cases, insulin-dependent diabetes patients generate antibodies against exogenous insulin. These antibodies can give rise to insulin-resistance and should therefore be controlled regularly when monitoring insulin therapy.

— The wells of the Anti-Insulin ELISA from ORGENTEC are coated with a mixture of bovine, pig, and recombinant human insulin.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 520	ORG 220-24
RANGE	0–100 U/ml	0–100 U/ml
CUT-OFF	10 U/ml	10 U/ml
CALIBRATION	internal	internal
COATING	bovine, pig, and recombinant human insulin	

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— beta-2-Microglobulin is an ELISA test for the quantitative determination of beta-2-microglobulin in human serum, plasma, or urine.

— beta-2-microglobulin, the light chain of the major histocompatibility complex class I (HLA I, MHC I) serves as a tumor marker in cases of multiple myeloma, Hodgkin's disease, chronic lymphatic leukemia, and other non-Hodgkin's lymphomas. Other diseases (infections, autoimmune diseases) that involve activation of the cellular immune system also cause increased beta-2-microglobulin concentrations in serum.

— beta-2-microglobulin is also used as a marker for evaluating the glomerular filtration rate, for the diagnosis and monitoring of tubule-interstitial kidney damage, for monitoring kidney function after a kidney transplant, and monitoring HIV infections.

— The wells of the microplate are coated with highly specific antibodies against human beta-2-microglobulin.

<b>SYSTEM</b>	<b>ELISA</b>
<b>TESTS</b>	<b>96 wells</b>
<b>CODE</b>	<b>ORG 5BM</b>
RANGE	0–12 µg/ml
CUT-OFF	3 µg/ml serum, 0.3 µg/ml urine
CALIBRATION	WHO B2M
COATING	anti-beta-2-microglobulin antibodies (polyclonal, rabbit)

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- Ferritin is an ELISA test for the quantitative determination of ferritin in human serum or plasma.
- It is used to diagnose iron deficiency and to monitor the treatment of high-risk individuals like pregnant women, blood donors, and hemodialysis patients.
- A further use for the determination of ferritin levels is in cases of iron overload (ferritin > 500 ng/ml), in individuals with conditions such as hemochromatosis and secondary siderosis. Extremely elevated ferritin levels are also found in cases of liver disease, neoplasia, infection, inflammation, and malignoma.
- The wells of the microplate in the Ferritin test are coated with highly specific antibodies against human ferritin.

<b>SYSTEM</b>	<b>ELISA</b>
<b>TESTS</b>	<b>96 wells</b>
<b>CODE</b>	<b>ORG 5FE</b>
RANGE	0–1500 ng/ml
CALIBRATION	WHO 94/572
COATING	anti-human-ferritin (polyclonal, rabbit)

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— Micro-Albumin is an ELISA test for the quantitative determination of human albumin in urine.

— The rate of albumin excretion is an important control and prevention parameter for the evaluation of kidney function in patients with diabetes. The Micro-Albumin ELISA test from ORGENTEC can be used to precisely determine even trace concentrations of albumin directly in urine.

— Normally, the concentrations of albumin in urine are very low. Values above 25 µg albumin per ml urine are an important indicator of kidney disease in diabetics. If the glomerular filtration process is highly active, albumin excretion may be elevated without any kidney damage. Physical stress may also cause increased albumin secretion.

SYSTEM	ELISA
TESTS	96 wells
CODE	ORG 5MA
RANGE	1.5–400 µg/ml
CUT-OFF	25 µg/ml
CALIBRATION	internal
COATING	albumin

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— *Bordetella pertussis*, the pathogen that causes whooping cough, is a small, Gram-negative bacterium in the coccobacillus family. It is transmitted by droplet infection and is very contagious.

— Whooping cough is one of the most dangerous infectious diseases for infants and can be fatal to them. In adolescents and adults the infection usually proceeds mildly, and is often not recognized as whooping cough. They can thus transmit the disease to vulnerable younger children.

— Laboratory diagnosis in infants can be accomplished through direct detection after bacterial culture or by PCR. Adolescents and adults mostly consult a doctor if symptoms persist for a longer period. Then, indirect serological methods are more suitable, such as an ELISA for the detection of antibodies against pertussis toxin, the critical pathogenic factor.

— The Anti-B. pertussis Toxin ELISA tests from ORGENTEC Diagnostika use highly purified native pertussis toxin as antigen for the detection of IgA or IgG antibodies. These two pure pertussis toxin ELISA tests allow for quantification. Both assays from ORGENTEC are calibrated according to the international WHO reference preparations.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 916G	ORG 916A
RANGE	0–400 IU/ml	0–60 IU/ml
CUT-OFF	50 IU/ml	12 IU/ml
CALIBRATION	WHO (NIBSC code 06/140 and 06/142)	WHO (NIBSC code 06/140 and 06/142)
COATING	purified native pertussis toxin	

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— The Anti-Borrelia tests for Alegria® allow for detection of IgG or IgM class antibodies against *Borrelia burgdorferi* sensu lato (s.l.), the main causative organisms of Lyme disease in Asia, Europe, and the USA. The tests use different sets of recombinant proteins displaying highly specific immunodominant epitopes.

— Both kits detect antibodies against VlsE (VMP-like sequence expression site) and OspC (outer-surface-protein C). In addition, Anti-Borrelia IgG (ORG 911G) detects IgG class antibodies against the chromosomally encoded high molecular weight protein p83/p100 and DbpA (decorin-binding-protein A). Anti-Borrelia IgM Abs. (ORG 911MX) contains the supplemental antigen p41i, the 14 kD internal flagellin fragment. These tailored antigen profiles cover the relevant epitopes for specific antibody recognition, avoiding cross-reactivity with other bacteria.

— The Anti-Borrelia IgM Abs. test is equipped with the ORGENTEC rheumatoid factor absorbent, which absorbs any rheumatoid factors that may be present during sample processing in Alegria®. This makes separate preparation of the serum sample unnecessary.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 911G	ORG 911MX
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	VlsE, DbpA, OspC, p83/p100	VlsE, OspC, p41i

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— The Anti-Borrelia Liquor tests for Alegria® are ELISA based test systems intended for the comparative quantitative measurement of IgG or IgM class antibodies against *Borrelia burgdorferi* sensu lato in human serum or plasma samples and in cerebrospinal fluid (CSF). The assays support the detection of specific antibody synthesis in the central nervous system.

— Neuroborreliosis is a disease that may develop when the bacteria cross the blood-brain barrier and infect the nervous system. A comprehensive CSF profil includes calculation of the specific anti-Borrelia antibody index (AI) as an evidence for *Borrelia* specific intrathecal antibody production. Calculation of AI is based on the determination of anti-Borrelia antibodies in a CSF sample and a corresponding serum sample from the same patient. The detection of an elevated *Borrelia* antibody index remains the main confirmatory tool in neuroborreliosis diagnostics.

— In these assays, recombinant highly specific immunodominant antigens of the three major human pathogenic *Borrelia* species are used as coating antigens. This leads to improved sensitivity and specificity of the ORGENTEC Anti-Borrelia Liquor ELISAs compared to tests based on purified antigens from cultivated bacteria.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 911GL	ORG 911ML
RANGE	10–200 U/ml	10–200 U/ml
CALIBRATION	internal	internal
COATING	VlsE, DbpA, OspC, p83/p100 (recombinant)	VlsE, OspC, p41i (recombinant)

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— The Anti-Chlamydia pneumoniae ELISA tests from ORGENTEC Diagnostika allow for the specific and automated detection of human antibodies against *Chlamydia pneumoniae*. These quantitative test systems support diagnosis of acute, latent, and chronic infections with this pathogen. The sample material is human serum or plasma.

— Diagnosis of *Chlamydia pneumoniae* infections is generally made after detection of specific antibodies in the serum.

— All ORGENTEC Anti-Chlamydia pneumoniae ELISAs are coated with elementary bodies and reticulate bodies of *Chlamydia pneumoniae*. By particular kind of antigen preparation, species specificity of these tests toward other species of *Chlamydia* is strikingly increased.

— The Anti-Chlamydia pneumoniae IgM Abs. test is equipped with the ORGENTEC rheumatoid factor absorbent, which absorbs any rheumatoid factors that may be present during sample processing in Alegria®. This makes separate preparation of the serum sample unnecessary.

SYSTEM	Alegria		
TESTS	24 Alegria® Test Strips		
CODE	ORG 907G	ORG 907A	ORG 907MX
RANGE	0–200 U/ml	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml	25 U/ml
CALIBRATION	internal	internal	internal
COATING	elementary and reticulate bodies of <i>Chlamydia pneumoniae</i> (CWL-029)		

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— The Anti-Chlamydia trachomatis ELISAs from ORGENTEC Diagnostika are test systems for the specific and automated detection of human IgA, IgG, and IgM antibodies against *Chlamydia trachomatis*. They can be used for the detection of acute, latent, or chronic infections with this bacterial pathogen and for determination of the immune status of the patient. The sample material is human serum or plasma.

— All ORGENTEC Anti-Chlamydia trachomatis ELISAs are coated with recombinant antigens of the major outer membrane protein complex (MOMP) of *Chlamydia trachomatis*. In addition two more recombinant antigens are used, TARP and CPAF, which bind *Chlamydia trachomatis*-specific antibodies. Thus, both test sensitivity and test specificity is raised compared to tests simply using MOMP antigens.

— The Anti-Chlamydia trachomatis IgM Abs. test is equipped with the ORGENTEC rheumatoid factor absorbent, which absorbs any rheumatoid factors that may be present during sample processing in Alegria®. This makes separate preparation of the serum sample unnecessary.

SYSTEM	Alegria		
TESTS	24 Alegria® Test Strips		
CODE	ORG 906G	ORG 906A	ORG 906MX
RANGE	0–200 U/ml	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml	25 U/ml
CALIBRATION	internal	internal	internal
COATING	CPAF, MOMP, TARP (all recombinant)		

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— The Anti-EBV (EBNA-1) IgG ELISA is used for the quantitative detection of human IgG antibodies against EBNA-1 of the Epstein-Barr virus. This ELISA is used for the detection of a latent infection. Samples used are serum or plasma.

— In the primary phase of infection, the pathogen first infects the salivary glands, which leads to general flu-like symptoms. After it moves to the B-lymphocytes, the virus spreads throughout the body, which can lead to the symptoms of infectious mononucleosis.

— About three weeks after the onset of symptoms, IgG antibodies against EBNA-1 are formed. Appearance of these antibodies indicates a shift from the active phase of the virus into a latent phase. IgG antibodies against EBNA-1 remain in the body throughout the life of the patient; however, reactivation of the virus can cause their concentration to drop below the detection limit of an ELISA.

— Acute EBV infections lead to the unspecific stimulation of the immune system, which can lead to false positive results in tests for other pathogens. It is thus always necessary to determine the patient's EBV status.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 902</b>
RANGE	0–800 U/ml
CUT-OFF	25 U/ml
CALIBRATION	internal
COATING	EBV nuclear antigen 1 (EBNA-1) (recombinant)

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— The ORGENTEC Anti-EBV (VCA) tests are quantitative ELISAs for the detection of specific antibodies against the Epstein-Barr virus in human serum. The Anti-EBV (VCA) IgM Abs. detects antibodies against the viral capsid antigen (VCA), which allows for the detection of an acute infection of the Epstein-Barr virus. The Anti-EBV (VCA) IgG is used to detect an existing acute or latent infection.

— In the early stage of the disease, IgM and IgG antibodies against VCA, the viral capsid antigen of the Epstein-Barr virus are detectable. About three weeks after the onset of symptoms, the maximum concentration of IgM antibodies against the VCA peptide is attained; the highest concentration of VCA IgG antibodies is reached after about six weeks. This high concentration of VCA IgG antibodies remains throughout the life of the patient.

— Anti-EBV (VCA) IgM Abs. is equipped with the ORGENTEC rheumatoid factor absorbent, which absorbs any rheumatoid factors that may be present during sample processing in Alegria®. This makes separate preparation of the serum sample unnecessary.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 901G	ORG 901MX
RANGE	0–400 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	EBV viral capsid peptides p18 and p23	

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— The Alegria® Anti-EBV (VCA) IgG Liquor is an ELISA based test system intended for the comparative quantitative measurement of IgG class antibodies against EBV (VCA) in human serum or plasma samples and in cerebrospinal fluid (CSF). This assay is used for the detection of IgG antibody synthesis in the central nervous system.

— The detection of pathogen-specific antibody production in CSF plays a role in the diagnosis of acute neurological infections. The antibody index (AI) reflects the specific immune reaction of one species of antibody. Calculation of the AI is based on testing for pathogen-specific antibodies in a CSF sample and a corresponding serum sample from the patient to detect specific intrathecal antibody production.

— Complications of an Epstein-Barr infection can result in a variety of neurological inflammatory processes (e.g. EBV encephalitis, meningitis, peripheral facial nerve paresis, acute myelopathies, radiculitis), though only in very rare cases. The detection of intrathecal IgG production against EBV is very useful for confirming a diagnosis, as well as for the differential diagnosis of neurological inflammatory processes of different etiology.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 901GL</b>
RANGE	10–200 U/ml
CALIBRATION	internal
COATING	EBV viral capsid peptide (VCAp18) (recombinant) , EBV viral capsid peptide (VCAp23) (recombinant)

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— Anti-EBV (ZEBRA) IgM is used for the quantitative determination of human IgM antibodies against certain peptide sequences in the ZEBRA protein (BZLF-1) of the Epstein-Barr virus. The detection of these IgM antibodies is used to support a diagnosis of acute primary infection with EBV. The sample material is human serum or plasma.

— It is generally possible to differentiate between the latent and lytic stages of infection with EBV. In the acute lytic phase of an infection all viral genes required for the formation of complete virus particles are activated. Included among the corresponding gene products is the ZEBRA protein (BamHI-Z-encoded Epstein Barr virus replication activator), an immediate early protein (IEA) of EBV. This is a DNA-binding nuclear protein that activates other genes; its expression precedes the synthesis of other molecules of the lytic phase.

— The target antigens of the Anti-EBV (ZEBRA) IgM test are IgM-specific peptides of the ZEBRA protein. The IgM antibodies it detects are found exclusively during the acute primary phase and are early markers of an EBV infection.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 910</b>
RANGE	0–200 U/ml
CUT-OFF	25 U/ml
CALIBRATION	internal
COATING	ZEBRA peptides

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— Anti-Helicobacter pylori IgA and Anti-Helicobacter pylori IgG are ELISA-based, automated test systems for the quantitative determination of IgA and IgG antibodies against *Helicobacter pylori* in human serum or plasma.

— About half of the global human population is infected with *H. pylori*. The infection is often acquired in early childhood and usually remains undetected. However, about 10–15 % of infected individuals develop severe stomach problems ranging from gastric ulcers to carcinoma of the stomach. Left untreated, *H. pylori* generally remains in the stomach throughout the life of the individual. Eradication is only possible by treatment with a combination of special antibiotics.

— There are a variety of invasive and non-invasive methods for detecting *H. pylori* infections. Non-invasive serological analyses are easy to carry out and have high sensitivity and specificity. Antibodies against the bacterial pathogenicity factors, particularly CagA are statistically more prevalent in cases of more serious infection. Specific antibodies against *H. pylori* are detectable by ELISA shortly after infection.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 917G	ORG 917A
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	<i>H. pylori</i> whole-cell lysate, recombinant CagA	

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— The hepatitis E virus is a small, non-enveloped virus with a single-stranded RNA genome. Along with hepatitis A, hepatitis E is one of the primary causes of acute hepatitis. The virus is transmitted by a faecal-oral route, through contaminated drinking water countries and through undercooked meat.

— Laboratory diagnosis of HEV infection depends on the detection of impaired liver function (elevated serum level of liver enzymes), as well as the direct or indirect detection of the virus. Because cases of hepatitis E cannot be distinguished clinically from other types of acute viral hepatitis, hepatitis E should always be considered in the differential diagnosis of acute hepatitis. During the viremic phase, HEV RNA can be detected in both the blood and stool by means of RT-PCR. However, the serological detection of virus-specific antibodies is more meaningful. The IgM response is generally positive four weeks after infection when clinical symptoms begin. The IgG antibodies then increase and reach their maximum about four weeks later.

— Anti HEV antibodies are predominantly directed against epitopes of ORF2. This test uses recombinant ORF2 from genotypes 1 and 3 as antigen. Due to cross reactions, antibodies against genotypes 2 and 4 are also detected.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 921G	ORG 921MX
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	recombinant ORF2 from genotypes 1 and 3	

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— The Anti-HSV-1 IgG and Anti HSV-1 IgM Abs. tests are ELISAs for the quantitative automated detection of human antibodies against type 1 herpes simplex viruses (HSV-1). Samples consist of serum or plasma. The Anti-HSV-1 IgM Abs. ELISA from ORGENTEC Diagnostika detects acute HSV-1 infections; Anti-HSV-1 IgG is used to determine the patient's immune status.

— The Anti-HSV-1 IgM Abs. is coated with purified HSV-1 antigen. This test enables the immediate and sensitive detection of a primary infection. The Anti-HSV-1 IgG ELISA use the type-specific glycoprotein G1, which allows for the differentiation of an HSV-1 infection from an infection with HSV-2. As a rule, IgG antibodies against HSV-1 are formed several weeks after a primary infection.

— IgM-class rheumatoid factors can distort the results of tests for the pathogen-specific detection of IgM antibodies. The Anti-HSV-1 IgM Abs. is equipped with the ORGENTEC rheumatoid factor absorbent that absorbs rheumatoid factors during sample processing in Alegria®. This makes prior serum sample preparation unnecessary.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 903G	ORG 903MX
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	HSV-1 glycoprotein G1	purified HSV-1 antigen

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— Anti-HSV-1/2 IgG and Anti-HSV-1/2 IgM Abs. by ORGENTEC Diagnostika are test systems for the detection of specific antibodies against human herpes simplex viruses (HSV). The Anti-HSV-1/2 IgM Abs. ELISA is designed to detect IgM antibodies against the viruses to identify acute infections. The detection of IgG antibodies against herpes simplex viruses is used to determine the immune status of patients. Test samples may be human serum or plasma.

— The Anti-HSV-1/2 assays from ORGENTEC Diagnostika are coated with antigen from inactivated type 1 and type 2 herpes simplex virus. The variety of epitopes allows for effective screening for IgG and IgM antibodies against both human herpes viruses. The type-specific Anti-HSV-1 and Anti-HSV-2 test systems from ORGENTEC Diagnostika are available for further differentiation of the antibody response.

— IgM-class rheumatoid factor antibodies can distort the results of tests for the pathogen-specific detection of IgM antibodies. The Anti-HSV-1/2 IgM Abs. (ORG 905MX) is equipped with a rheumatoid factor absorbent that absorbs rheumatoid factors during sample processing in Alegria®. This makes prior serum sample preparation unnecessary.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 905G	ORG 905MX
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	purified HSV-1 and HSV-2 antigens	

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— Anti-HSV-1/2 IgG Liquor is an ELISA based test system intended for the comparative quantitative measurement of IgG class antibodies against HSV-1/2 in human serum or plasma samples and in cerebrospinal fluid (CSF). This assay is used for the detection of IgG antibody synthesis in the central nervous system (CNS).

— The detection of pathogen-specific antibody production in CSF plays a role in the diagnosis of acute neurological HSV infections. The antibody index (AI) reflects the ratio of HSV-specific antibodies between CSF and serum. Calculation of the AI is based on the detection of pathogen-specific antibodies in a CSF sample and a corresponding serum sample from the patient to detect specific intrathecal antibody production.

— Inflammation of the CNS may occur as rare complication of infection with HSV, in particular HSV-associated encephalitis, but also meningitis, peripheral facial nerve palsy, myelitis, and radiculitis. In cases of HSV encephalitis, production of antibodies in the CSF begins in the second week and slowly decreases over months and years. Intrathecal production of IgG against HSV is also often detected in chronic inflammatory processes like multiple sclerosis.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 905GL</b>
RANGE	10–200 U/ml
CALIBRATION	internal
COATING	HSV-1 whole-virus antigen, HSV-2 whole-virus antigen

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— The Anti-HSV-2 IgG and Anti-HSV-2 IgM Abs. test systems detect human antibodies against type 2 herpes simplex viruses (HSV-2). Anti-HSV-2 IgM Abs. is for the automated detection of acute HSV-2 infection. Anti-HSV-2 IgG ELISA is used to determine the immune status of the patient.

— For the immediate and sensitive detection of a primary infection, the Anti-HSV-2 IgM Abs. ELISA is coated with purified HSV-2 antigen. For the detection of IgG antibodies, the Anti-HSV-2 IgG test uses the type-specific glycoprotein G2.

— As a rule, IgG-class antibodies against glycoprotein G2 are detectable several weeks after primary infection. Their detection makes it possible to differentiate between an infection with HSV-1 and HSV-2.

— IgM-class rheumatoid factor antibodies can distort the results of tests for the pathogen-specific detection of IgM antibodies. The Anti-HSV-2 IgM Abs. is equipped with the ORGENTEC rheumatoid factor absorbent that absorbs rheumatoid factors during sample processing in Alegria®. This makes prior serum sample preparation unnecessary.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 904G	ORG 904MX
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	HSV-2 glycoprotein G2	purified HSV-2 antigens

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— Anti-Measles Virus IgG and Anti-Measles Virus IgM Abs. are ELISA-based test systems for the quantitative determination of IgG or IgM antibodies against measles virus in human serum or plasma.

— IgM antibodies against the measles virus can be detected about two weeks after primary infection for a period of four to six weeks. They can be detected and quantified with the Anti-Measles-Virus IgM assay. IgM-class rheumatoid factor antibodies can distort the results of pathogen-specific IgM antibody tests. The Anti-Measles Virus IgM test is equipped with the ORGENTEC rheumatoid factor absorbent. This reagent absorbs any rheumatoid factors that may be present during sample processing in the Alegria® instrument, making separate preparation of the serum sample unnecessary.

— Once infected, a patient generally develops life-long immunity, measles virus IgG antibodies remain detectable throughout the lifetime of the patient and indicate the patient's immune status. Calibration of the Anti-Measles Virus IgG test by using the WHO NIBSC 97/648 standard facilitates comparison of results across different laboratories.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 909G	ORG 909MX
RANGE	0–5000 mIU/ml	0–200 U/ml
CUT-OFF	250 mIU/ml	25 U/ml
CALIBRATION	WHO 97/648	internal
COATING	purified particles of inactivated measles virus (strain Edmonston)	

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— Anti-Measles Virus IgG Liquor is an ELISA-based, automated system for the comparative quantitative determination of IgG antibodies against measles virus in human serum or plasma and in cerebrospinal fluid (CSF). This test is used to detect IgG antibody synthesis in the central nervous system.

— Measles infection causes encephalitis in about one in 1000 children. Subacute sclerosing panencephalitis (SSPE) occurs in approximately one in 10000 children and is always fatal. An increased titer of specific antibodies to measles virus is found in the CSF of the affected children.

— In chronic diseases with CNS involvement, e.g. multiple sclerosis (MS), an elevated antibody index of anti-measles virus, anti-rubella virus and/or anti-varicella zoster virus (MRZ reaction) can be detected in the CSF. In about 90 % of MS patients, an increased antibody index against the MRZ combination can be detected.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 909GL</b>
RANGE	10–200 U/ml
CALIBRATION	internal
COATING	purified particles of inactivated measles virus (strain Edmonston)

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— Anti-Mumps Virus IgG and IgM Abs. by ORGENTEC Diagnostika are quantitative ELISA test systems for the specific detection of IgG or IgM antibodies against mumps virus in human serum or plasma. The IgM Abs. ELISA is used for the detection of acute infections; the IgG ELISA is used to determine the immune status of a patient.

— In most patients, IgM antibodies against mumps virus can be detected shortly after the appearance of clinical symptoms; sometimes they can persist for several months.

— Mumps viruses are present around the world. Humans are their only natural host. Once an infection, either clinically manifested or symptomless, has been overcome, the patient mostly develops life-long immunity.

— IgM class rheumatoid factors may lead to false results in tests for the pathogen-specific detection of IgM antibodies. The Anti-Mumps Virus IgM Abs. features the ORGENTEC rheumatoid factor absorbent that absorbs rheumatoid factors during sample processing in Alegria®. No preparation of the serum sample is thus required.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 908G	ORG 908MX
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	purified particles of inactivated mumps virus (strain Enders)	

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— *Mycoplasma pneumoniae* is one of the smallest free-living organisms. It forms no cell wall and has limited biosynthetic capabilities. In addition to *Chlamydia* and *Legionella*, *M. pneumoniae* is one of the primary causes of community-acquired pneumonia. Differentiation of the pathogens is crucial for choosing the suitable antibiotics.

— Serological detection of pathogen-specific antibodies with ELISA is a common method used to diagnose infections caused by *M. pneumoniae*. Antibodies against *M. pneumoniae*'s adhesins are particularly useful for this. These antibodies are characterized by their high sensitivity and specificity.

— The individual detection of IgG, IgM, and IgA antibodies is sufficient to differentiate between primary infections and reinfections. IgM antibodies appear seven to ten days after onset of the infection, they are often absent in reinfections. A significant increase in the IgG titre within two weeks is indicative of an acute infection. High IgA titres are the best indicator of acute infection in adults; however low titres of IgA antibodies can persist over longer periods of time.

SYSTEM	Alegria		
TESTS	24 Alegria® Test Strips		
CODE	ORG 915G	ORG 915A	ORG 915MX
RANGE	0–200 U/ml	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml	25 U/ml
CALIBRATION	internal	internal	internal
COATING	<i>M. pneumoniae</i> whole-cell lysate, adhesin P1 (recombinant)		

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— Parvovirus B19 forms two capsid proteins VP1 and VP2, as well as the non-structural protein NS1. In the acute phase of the infection IgG and IgM antibodies bind to linear epitopes of the entire structural protein and to conformational epitopes of the virus capsid. After decline of the initial IgM and IgG response IgG antibodies persist, which are preferentially directed against conformational epitopes.

— Recombinant VP2 virus-like particles are the antigens in both tests. Presentation of the correct tertiary structure of the antigens is prerequisite for the detection of IgG antibodies in the acute phase of the infection and in past infections.

— Anti-Parvovirus B19 IgM Abs. detects IgM class antibodies, allowing for identification of a recent primary infection. The Parvovirus B19 IgM Abs. test is equipped with the ORGENTEC rheumatoid factor absorbent, which absorbs any rheumatoid factors that may be present during sample processing in Alegria®. Separate preparation of the serum sample is not necessary.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 912G	ORG 912MX
RANGE	0–50 IU/ml	0–200 U/ml
CUT-OFF	3 IU/ml	25 U/ml
CALIBRATION	WHO (NIBSC code 01/602)	internal
COATING	VP2 virus-like particles (recombinant)	

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— Anti-Rubella Virus IgG Liquor is an ELISA based test system intended for the comparative quantitative measurement of IgG class antibodies against rubella Virus in human serum or plasma samples and in cerebrospinal fluid (CSF). This assay is applicable for the detection of IgG antibody synthesis in the central nervous system.

— Rubella infection sometimes leads to serious neurological complications and causes severe brain injury in children if their mother is infected in early pregnancy. Rubella virus can also cause a progressive “pan-encephalitis” later in life in children who were infected in the womb and survived but remain chronically infected.

— In chronic diseases with CNS involvement, e.g. multiple sclerosis (MS) an elevated antibody index of anti-measles virus, anti-rubella virus and/or anti-varicella zoster virus can be detected in the CSF (MRZ reaction). In about 90 % of MS patients, an increased antibody index against the MRZ combination can be detected.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 919GL</b>
RANGE	10–200 U/ml
CALIBRATION	internal
COATING	purified particles of inactivated rubella virus (strain HPV-77)

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— Varicella zoster virus (VZV), cause of chicken pox and shingles (Herpes zoster), belongs to the family of alpha herpes viruses that are pathogenic to humans. Transmission occurs by droplet and smear infection or transplacental. After an incubation period of two weeks fever and the typical itchy chicken pox rash occur. Primary infection during the first 20 weeks of gestation, may cause severe abnormalities in the foetus. Neonatal chicken pox infection may also be quite severe.

— Chicken pox and Herpes zoster are usually diagnosed by means of their clinical symptoms. In cases with complications or infections during pregnancy, laboratory diagnosis is important. IgM antibodies are typical of an acute infection. Reactivation of the virus (Herpes zoster) induces an increase of IgA and IgG antibody titres.

— The tests are based on VZV glycoproteins as antigens. The Anti-VZV Abs. test is equipped with the ORGENTEC rheumatoid factor absorbent, which absorbs any rheumatoid factors that may be present during sample processing in Alegria®. Separate preparation of the serum sample is not necessary.

SYSTEM	Alegria		
TESTS	24 Alegria® Test Strips		
CODE	ORG 914G	ORG 914A	ORG 914MX
RANGE	0–1500 mIU/ml	0–200 U/ml	0–200 U/ml
CUT-OFF	100 mIU/ml	25 U/ml	25 U/ml
CALIBRATION	WHO (NIBSC code W1044)	internal	internal
COATING	purified VZV glycoproteins from virus cell culture		VZV glycoprotein E (recombinant)

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— Anti-VZV IgG Liquor is an ELISA-based, automated test system for the comparative quantitative determination of IgG antibodies against varicella zoster virus (VZV) in human serum or plasma and in cerebrospinal fluid (CSF). It is used for the detection of IgG antibody synthesis in the central nervous system.

— An elevated varicella zoster virus (VZV) antibody index characterizes varicella zoster ganglionitis, often accompanied by facial paralysis. VZV-antibody detection is important for the differential diagnosis of neuroborreliosis. CSF analysis currently is the method of choice to distinguish these infections that may both be the cause of facial paralysis.

— In chronic diseases with CNS involvement, e.g. multiple sclerosis (MS) an elevated antibody index of anti-measles virus, anti-rubella virus and/or anti-varicella zoster virus can be detected in the CSF (MRZ reaction). In about 90 % of MS patients, an increased antibody index against the MRZ combination can be detected.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 914GL</b>
RANGE	10–200 U/ml
CALIBRATION	internal
COATING	VZV glycoproteins

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— *Yersinia* are found worldwide in animal reservoirs. Infections with *Y. enterocolitica* and *Y. pseudotuberculosis* are common. *Y. enterocolitica* initially infects the gastrointestinal tract, causing fever, vomiting, and bloody diarrhoea. After two to four weeks reactive arthritis may develop. Infections with *Y. pseudotuberculosis* may lead to a condition known as pseudoappendicitis, which is often mistaken for true appendicitis.

— Because *Yersinia* are only detectable in stool for a short time and are not detectable at all in the joints affected by arthritis, the serological detection of *Yersinia*-specific antibodies in serum is important. The ORGENTEC test uses recombinant YopD for the detection of antibodies against *Yersinia* sp. The Yop proteins (*Yersinia* outer membrane proteins) are plasmid-encoded virulence factors that are highly specific for all of the *Yersinia* species.

— After initial contact with this pathogen, IgM antibodies are formed, and disappear again shortly after. IgG and IgA antibodies follow. High IgA and IgG antibody titres are indicative of a persistent infection as a possible cause of arthritis.

SYSTEM	Alegria	
TESTS	24 Alegria® Test Strips	
CODE	ORG 913G	ORG 913A
RANGE	0–200 U/ml	0–200 U/ml
CUT-OFF	25 U/ml	25 U/ml
CALIBRATION	internal	internal
COATING	recombinant YopD	

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— Alpha-1-Antitrypsin for Alegria® is an ELISA-based, automated, in-vitro test system for the quantitative determination of Alpha-1-Antitrypsin in stool. It is used to evaluate the activity level of chronic inflammatory bowel diseases. In addition, faecal alpha-1-antitrypsin is also a marker for protein leakage or increased permeability of the intestinal mucosa when the mucosal barrier is not intact.

— This protein is a primary inhibitor that reversibly forms complexes with serine proteases in the digestive tract and with active inflammatory immune cells. Alpha-1-antitrypsin plays an important regulatory and anti-inflammatory role, preventing an attack on neighbouring healthy tissue during an inflammatory event. This function makes alpha-1-antitrypsin an indicator for active inflammatory processes.

— Elevated values for alpha-1-antitrypsin in stool are primarily found in cases of ulcerative colitis and Crohn's disease. However, they are also found in cases of enteral protein leakage from other causes, such as Whipple's disease, necrotising enterocolitis, intestinal tuberculosis, celiac disease, intestinal lymphoma, or systemic lupus erythematosus. The levels of alpha-1-antitrypsin may also be elevated in patients with intestinal polyps, colon cancer, diverticulitis, food allergies, and in pregnancy.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 286</b>
RANGE	80–3600 µg/g
CALIBRATION	alpha-1-antitrypsin reference preparation
COATING	anti-alpha-1-antitrypsin antibodies

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— Calprotectin is an ELISA test for the quantitative determination of calprotectin in stool samples. The test contributes to the evaluation of inflammatory bowel disease (IBD) and to the differentiation of IBD and irritable bowel syndrome (IBS).

— Non-specific inflammatory parameters in blood and clinical symptoms are not sufficiently specific to diagnose inflammation localized in the intestines. Diagnosis of inflammatory bowel diseases is thus often made with considerable delay. The detection of calprotectin in stool samples complements basic laboratory diagnostic tests with a highly diagnostic intestinal inflammatory parameter. Elevated calprotectin levels in stool are found in cases of chronic inflammatory bowel disease, infectious gastroenteritis and colorectal carcinoma, among others. In patients with irritable bowel syndrome, calprotectin levels are in the normal range.

— Determination of fecal calprotectin is a simple, non-invasive test that can easily be carried out before any imaging or endoscopic procedures, and aids in monitoring progression of the disease and its treatment.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 580	ORG 280
RANGE	15-1000 µg/g	0-1000 µg/g
CUT-OFF	50 µg/g	50 µg/g
CALIBRATION	internal	internal
COATING	anti-calprotectin antibodies	

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— Lactoferrin is an ELISA-based test system for the quantitative measurement of Lactoferrin in human stool samples, to be used in the assessment of intestinal inflammatory disorders.

— Non-specific inflammatory parameters in blood (CRP, ESR, complete blood count) combined with clinical symptoms are not sufficiently specific to diagnose inflammation localized in the intestines. The detection of lactoferrin in stool samples complements basic laboratory diagnostic tests with a highly significant intestinal inflammatory parameter.

— Lactoferrin is an iron-binding protein present in polymorphonuclear granulocytes. The fecal Lactoferrin level correlates with the intensity of the inflammation of the intestinal mucosa.

— Lactoferrin in stool is a sensitive marker for inflammatory bowel disease (IBD) and aids in monitoring progression of the disease and its treatment.

<b>SYSTEM</b>	<b>Alegria</b>
<b>TESTS</b>	<b>24 Alegria® Test Strips</b>
<b>CODE</b>	<b>ORG 284</b>
RANGE	0–150 µg/g
CUT-OFF	7.2 µg/g
CALIBRATION	internal
COATING	anti-lactoferrin antibodies

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— 25-OH Vitamin D<sub>3</sub>/D<sub>2</sub> is an ELISA-based, automated, in-vitro test system for the quantitative determination of the total concentration of 25-OH vitamin D<sub>3</sub> and 25-OH vitamin D<sub>2</sub> in human serum or plasma (EDTA plasma, heparin plasma, citrate plasma) for the evaluation of vitamin D levels. This test uses an antibody that binds both forms of vitamin D equally well and with high affinity.

— Vitamin D<sub>3</sub> (cholecalciferol) is produced in the skin under irradiation with UV-B rays from the sun. It is metabolized to 25-OH vitamin D<sub>3</sub> in the liver. Vitamin D<sub>2</sub> (ergocalciferol) is mostly produced synthetically and is used as a dietary supplement. It is also converted to its corresponding 25-OH derivative by the liver.

— For reliable determination of the total vitamin D concentration, a special sample buffer is used to release the 25-OH vitamin D<sub>3</sub> and 25-OH vitamin D<sub>2</sub> from the vitamin D binding protein complex in the serum. No separate sample preparation steps are required.

SYSTEM	ELISA	Alegria
TESTS	96 wells	24 Alegria® Test Strips
CODE	ORG 570	ORG 270
RANGE	5–120 ng/ml	5–200 ng/ml
CUT-OFF	20 ng/ml (50 nmol/l)	20 ng/ml (50 nmol/l)
CALIBRATION	REF!25OHD	REF!25OHD
COATING	anti-25-OH Vitamin D <sub>3</sub> /D <sub>2</sub> antibodies	

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— 25-OH Vitamin D<sub>3</sub>/D<sub>2</sub> Control Set contains external control materials for the quantitative Alegria® 25-OH Vitamin D<sub>3</sub>/D<sub>2</sub> assay. The control set includes 2.5 ml each of Control A and Control B with a defined concentration of 25-OH vitamin D<sub>3</sub>/D<sub>2</sub>. Set is for 30 determinations

— These controls are available separately.

<b>SYSTEM</b>	<b>Control Set</b>
<b>TESTS</b>	<b>2 x 2.5 ml</b>
<b>CODE</b>	<b>ORG 271</b>

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— The ANA-9-Line immunoblot test is a membrane-based enzyme immunoassay for the semi-quantitative determination of IgG antibodies against nine different nuclear antigens (ANA) in human serum or plasma. This immunoblot test is used for the differential diagnosis of rheumatic diseases.

— The nine purified antigens are deposited as bands onto nitrocellulose strips. The bound antibodies are stained by an enzyme-catalysed substrate reaction for detection.

— The ANA-9-Line blot is an important and easy-to-use test system for the detection of autoantibodies typical of various rheumatic autoimmune diseases; some of these autoantibodies act as markers for individual diseases. This test is particularly valuable for the differentiation of the autoantibody profile after ELISA screening or immunofluorescence tests.

SYSTEM	IMMUNOBLOT	
	8 strips	16 strips
TESTS	ORG 710-08	ORG 710-16
CODE		
EVALUATION	semiquantitative	semiquantitative
CUT-OFF	calibration strip	calibration strip
CALIBRATION	internal	internal
COATING	SS-A 52 (Ro 52), SS-A 60 (Ro 60), SS-B (La), Sm, RNP/Sm, Scl-70, Jo-1, centromere B, Rib-P	

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— Myositis plus is a membrane-based enzyme immunoassay for the semi-quantitative determination of IgG autoantibodies in human serum or plasma. This immunoblot test is used in cases of suspected polymyositis/dermatomyositis or overlap syndromes as well as for the differential diagnosis of scleroderma.

— Autoantibodies against nuclear and/or cytoplasmic antigens are found in nearly 90 % of patients with poly- or dermatomyositis. The Myositis plus test detects myositis-specific autoantibodies that are generally only found in cases of idiopathic inflammatory myopathy and thus act as markers for this condition. In addition, this test detects myositis-associated autoantibodies, which are detected in patients with other diseases as well as those with myopathy.

— The nine antigens in the Myositis plus test are deposited in bands onto nitrocellulose strips. The test detects autoantibodies characteristic of autoimmune myositis, scleroderma, and overlap syndromes in a single and simple step.

SYSTEM	IMMUNOBLOT	
	8 strips	16 strips
TESTS	8 strips	16 strips
CODE	ORG 760-08	ORG 760-16
EVALUATION	semiquantitative	semiquantitative
CUT-OFF	calibration strip	calibration strip
CALIBRATION	internal	internal
COATING	M2 proteins (PDC-E2, BCOADC-E2, OGDC-E2), Jo-1, PM/Scl-100, PL-7, PL-12, Mi-2, Ku (p70/80), SRP, Rib-P	

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- The Nucleo-9-Line test is a membrane-based enzyme immunoassay for the semi-quantitative determination of IgG autoantibodies against various nuclear antigens (ANA) in human serum or plasma.
- The immunoblot Nucleo-9-Line uses nine antigens that typically appear in various inflammatory rheumatologic conditions. Especially Sm, dsDNA and nucleosomes play an important role in the diagnosis of systemic lupus erythematosus (SLE). The nine purified antigens are individually applied as bands on a nitrocellulose strip. For detection bound antibodies are stained by an enzyme-catalysed substrate reaction.
- The Nucleo-9-Line blot is an important and easy-to-use test system for the further differentiation of ANA after ELISA screening or immunofluorescence tests.

SYSTEM	IMMUNOBLOT	
	8 strips	16 strips
TESTS	ORG 711-08	ORG 711-16
CODE	ORG 711-08	ORG 711-16
EVALUATION	semiquantitative	semiquantitative
CUT-OFF	calibration strip	calibration strip
CALIBRATION	CDC sera	CDC sera
COATING	dsDNA, nucleosomes, SS-A (Ro 52, Ro 60), SS-B (La), Sm, RNP/Sm, Scl-70, Jo-1, centromere B	

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— The ANCA-3-Line immunoblot test is a membrane-based enzyme immunoassay for the semi-quantitative determination of antineutrophil cytoplasmic IgG antibodies (ANCA) against proteinase 3 (PR3), myeloperoxidase (MPO), and glomerular basement membrane (GBM) in human serum or plasma. It is used for diagnostics of systemic autoimmune vasculitis.

— ANCA associated vasculitis may be a fast progressing life threatening disease, thus rapid and reliable ANCA serology is particularly vital in vasculitis diagnosis. ANCA-3-Line is a robust and efficient test system for rapid simultaneous analysis of characteristic antibody specificities for granulomatosis with polyangiitis (GPA, Wegener's granulomatosis), microscopic polyangiitis (MPA) and Goodpasture syndrome.

— Verification of the blot results and subsequent quantitative determination of the three individual ANCA can be carried out with specific ELISA tests: Anti PR3 *hs* (*high sensitive*), Anti-MPO, and Anti-GBM.

SYSTEM	IMMUNOBLOT	
	8 strips	16 strips
TESTS	8 strips	16 strips
CODE	ORG 789-08	ORG 789-16
EVALUATION	semiquantitative	semiquantitative
CUT-OFF	calibration strip	calibration strip
CALIBRATION	internal	internal
COATING	proteinase 3 (PR3), myeloperoxidase (MPO), glomerular basement membrane (GBM)	

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- The Gastro-5-Line immunoblot test is a membrane-based enzyme immunoassay for the semi-quantitative determination of IgG and IgA autoantibodies in human serum or plasma when autoimmune disease of the gastrointestinal tract is suspected. Gastro-5-Line used as an aid in the differential diagnosis of inflammatory autoimmune diseases of the digestive tract.
- Each nitrocellulose strip is coated with the antigens most commonly found in cases of gastrointestinal autoimmune disease:
  - Antibodies against mannan from *Saccharomyces cerevisiae* (ASCA) are ideal for the differentiation of Crohn's disease and ulcerative colitis.
  - Anti-tissue transglutaminase is a reliable marker for celiac disease. In addition, celiac patients also have serum antibodies against gliadin. Celiac patients often form no IgA antibodies. For this reason, the Gastro-5-Line immunoblot simultaneously detects both IgG and IgA class antibodies.
  - Antibodies against intrinsic factor are found in patients with pernicious anemia. Because of its high sensitivity, parietal cell antigen is a good screening parameter for pernicious anemia.

SYSTEM	IMMUNOBLOT	
TESTS	8 strips	16 strips
CODE	ORG 740-08	ORG 740-16
EVALUATION	semiquantitative	semiquantitative
CUT-OFF	calibration strip	calibration strip
CALIBRATION	internal	internal
COATING	mannan ( <i>Saccharomyces cerevisiae</i> ), tissue transglutaminase (tTG), gliadin, intrinsic factor, H <sup>+</sup> /K <sup>+</sup> -ATPase from parietal cells	

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— Liver-9-Line 2nd Generation is a membrane-based enzyme immunoassay for semi-quantitative determination of IgG autoantibodies against nine test antigens in human serum or plasma, to provide a differentiated serological profile for autoimmune liver diseases.

— Elevated titers of antinuclear antibodies (ANA) and/or anti-smooth muscle antibodies (ASMA) are typical of type 1 autoimmune hepatitis (AIH-1). Detection of antibodies against the individual smooth-muscle antigens F-Actin, alpha-actinin, and tropomyosin serves to differentiate this disease, for assessment of disease activity and evaluation of further progression. Autoantibodies against soluble liver antigen (SLA) are also highly specific markers for AIH-1. Antibodies against liver-kidney microsome type 1 (LKM-1) and liver cytosolic protein type 1 (LC1) are characteristic of AIH-2.

— Anti-mitochondrial antibodies of the M2 subtype (AMA-M2) indicate PBC. If they are absent the highly specific antibodies against nuclear body antigen Sp100 and nuclear pore antigen gp210 are diagnostically important.

SYSTEM	IMMUNOBLOT	
	8 strips	16 strips
TESTS	8 strips	16 strips
CODE	ORG 722-08	ORG 722-16
EVALUATION	semiquantitative	semiquantitative
CUT-OFF	calibration strip	calibration strip
CALIBRATION	internal	internal
COATING	F-actin, alpha-actinin, tropomyosin, SLA/LP, LKM-1, LC1, M2 proteins (PDC-E2, BCOADC-E2, OGDC-E2), Sp100, gp210	

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— **rheumachec® rapid test** is a membrane-based lateral flow assay for the rapid combined detection of rheumatoid factors (RF) and autoantibodies against mutated citrullinated vimentin (Anti-MCV®) in whole blood in cases of suspected rheumatoid arthritis.

— **rheumachec® rapid test** assists the doctor in the detection of these serological markers for early rheumatoid arthritis immediately at the point of care, at the first suspicion of rheumatoid arthritis. The test requires no laboratory equipment.

— A drop of capillary blood from the patient is merely applied to the test strip. A few drops of buffer solution start the test running, and after only 15 minutes, the test result shows if MCV® antibodies or rheumatoid factors are present in the blood sample. An abnormal test result should be followed up by immediate referral to a rheumatologist for further diagnostics and rapid treatment of the patient.

— In cases of rheumatoid arthritis, making treatment decisions early is particularly critical to successful therapy. When patients with rheumatoid arthritis are diagnosed and treated in time permanent damage can be avoided. Sometimes the disease can be put into long-term remission.

SYSTEM	Lateral flow test	
TESTS	1 Test	10 Tests
CODE	ORG 170-01	ORG 170-10
RANGE	qualitative	qualitative
CALIBRATION	internal	internal
COATING	MCV (mutated citrullinated vimentin), Fc fragments of IgG	

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— Tissue sections from monkey esophagus are the substrate for detection of Anti-Endomysium antibodies (AEA). AEA of the immunoglobulin A class (IgA) are found in patients with active celiac disease or dermatitis herpetiformis. Their diagnostic sensitivity for celiac disease amounts to 100 % in untreated patients, who do not adhere to a gluten free diet; the diagnostic specificity of AEA is 99 to 100 %.

— AEA may be used for disease monitoring. About twelve months after starting with the gluten free diet recommended for treatment of celiac disease, the levels of endomysial antibodies should drop below the detection limit.

— International guidelines for diagnosis of celiac disease approve testing for IgA class antibodies against tissue transglutaminase 2 (TG2, tTG) for initial diagnosis. Performing indirect immunofluorescence tests for AEA is always recommended if these first line tests give ambiguous results.

<b>SYSTEM</b>	<b>IFT KIT</b>
<b>TESTS</b>	<b>96 (12 slides, 8 wells each)</b>
<b>CODE</b>	<b>ORG 872</b>
COMPONENTS	Anti-Endomysium Antibodies (AEA) slides; AEA Positive Control IgA FITC/Evans; Mounting Medium PBS (10×); Negative Control; Blotting Paper 1×8

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— Detection of antibodies against native double-stranded DNA (Anti-nDNA, Anti-dsDNA) is very important for laboratory diagnostics of systemic lupus erythematosus (SLE). Immunofluorescence assays based on antigens of the flagellate *Crithidia luciliae* have a high diagnostic specificity. Sensitivity is somewhat lower, especially in patients with no signs for an active disease state.

— Antibodies against native DNA are more than specific markers for the disease; they are also suitable parameters for measuring lupus disease activity. Therefore, Anti-nDNA antibodies are among the most frequently detected autoantibodies that are associated with SLE. They are found in 95 % of SLE patients with renal involvement, in 50–70 % of SLE patients without renal involvement and in 40 % of patients with inactive SLE. Anti-nDNA antibodies are rarely found in healthy individuals.

SYSTEM	IFT KIT
TESTS	120 (10 slides, 12 wells each)
CODE	ORG 871
COMPONENTS	Anti-nDNA Antibodies (nDNA) slides; nDNA Positive Control; IgG FITC/Evans; Mounting Medium; PBS (10 ×); Negative Control; Blotting Paper 1 × 12

# Immunofluorescence

## Anti-Neutrophil Cytoplasmic Antibodies (ANCA-Ethanol)

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— Autoantibodies against cytoplasmic antigens of neutrophil granulocytes (ANCA) are acknowledged serological markers for the diagnosis of systemic vasculitis. Our slides contain ethanol-fixed neutrophil granulocytes, which facilitate the differentiation of granulomatosis with polyangiitis (GPA; formerly Wegener’s granulomatosis) from eosinophilic granulomatosis with polyangiitis (EGPA; formerly Churg-Strauss syndrome), microscopic polyangiitis (MPA), or Goodpasture syndrome.

— A granular immunofluorescence pattern in the cytosol (cANCA) indicates antibodies against proteinase 3 (PR3), which are indicative of GPA. A dense immunofluorescence pattern around the cell nucleus (pANCA) denotes antibodies against myeloperoxidase (MPO), which occur in cases of EGPA or MPA. Our formalin-fixed granulocyte test ORG 892 should be used subsequent to a pANCA positive test result with ethanol-fixed granulocytes to differentiate a true pANCA pattern from an ANA result with a similar perinuclear fluorescence pattern.

— In accordance with the “ANCA Consensus Statement” we recommend further serological analysis with an antigen-specific ELISA test for the individual detection of anti-PR3 (ORG 618, ORG 318-24) and anti-MPO (ORG 519, ORG 219-24).

SYSTEM	IFT KIT
TESTS	60 (10 slides, 6 wells each)
CODE	ORG 891
COMPONENTS	Anti-Neutrophil Cytoplasmic Antibodies (ANCA-Ethanol) slides; cANCA Positive Control; pANCA Positive Control; IgG FITC/Evans; Mounting Medium; PBS (10x); Negative Control; Blotting Paper 1x6

## Immunofluorescence

### Anti-Neutrophil Cytoplasmic Antibodies (ANCA-Formalin)

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— Indirect immunofluorescence tests for the detection of autoantibodies against cytoplasmic antigens of neutrophil granulocytes (ANCA) facilitate the serological diagnosis of autoimmune small-vessel vasculitis. Our formalin-fixed slides should be used subsequent to a pANCA-positive test result with ethanol-fixed granulocytes. They allow for the differentiation of a true pANCA pattern from antinuclear antibodies (ANA), which result in a similar perinuclear fluorescence pattern on ethanol-fixed granulocytes (ORG 891).

— If a granular fluorescence pattern is detected in the cytosol with our formalin-fixed slides, antibodies against myeloperoxidase are present. These may be found in cases of eosinophilic granulomatosis with polyangiitis (EGPA; formerly Churg-Strauss syndrome), microscopic polyangiitis (MPA), or Goodpasture syndrome.

— In accordance with the “ANCA Consensus Statement” we recommend further serological diagnostic analysis of the cytosolic (cANCA) or perinuclear (pANCA) immunofluorescence pattern with an antigen-specific ELISA test for the detection of anti-PR3 (ORG 618, ORG 318-24) and anti-MPO (ORG 519, ORG 219-24).

SYSTEM	IFT KIT
TESTS	60 (10 slides, 6 wells each)
CODE	ORG 892
COMPONENTS	Anti-Neutrophil Cytoplasmic Antibodies (ANCA-Formalin) slides; pANCA Positive Control; IgG FITC/Evans; Mounting Medium; PBS (10×); Negative Control; Blotting Paper 1×6

# Immunofluorescence

## Anti-Nuclear Antibodies HEp-2 (ANA-HEp-2)

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— Our cultured Hep-2 cells allow for simultaneous detection of a variety of autoantibodies with a single analytical preparation. Evaluation of different fluorescence patterns of the nucleus and the cytoplasm provides additional information about the bound antibody.

— Homogenous staining of the cell nucleus is indicative of systemic lupus erythematosus (SLE). A speckled pattern is strongly related to SLE, but it may also occur in patients with mixed connective tissue disease, Sjogren's syndrome, polymyositis or scleroderma. Peripheral fluorescence can be observed in patients with autoimmune liver diseases. A nucleolar pattern has been described in patients with overlapping scleroderma and polymyositis/dermatomyositis syndromes and with systemic scleroderma, especially for those with renal involvement. In patients with the limited cutaneous form of systemic sclerosis staining of the centromeres can be detected.

SYSTEM	IFT KIT
TESTS	120 (10 slides, 12 wells each)
CODE	ORG 870
COMPONENTS	Anti-Nuclear Antibodies HEp-2 (ANA-HEp-2) slides; ANA-HEp-2 Positive Control; IgG FITC/Evans; Mounting Medium; PBS (10x); Negative Control; Blotting Paper 1 x 12

# Immunofluorescence

## Anti-Smooth Muscle Antibodies (ASMA)

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— Tissue sections from rat stomach are used for the detection of anti-smooth muscle antibodies (ASMA) in serum or plasma. ASMA bind to the corresponding antigens in the muscle layers of the blood vessels and the interglandular fibrils in the rat stomach.

— Differentiation of different subgroups of ASMA antibodies, ASM-V, ASMA-G and ASMA-T, is made by evaluation of vessels, glomeruli and peritubular fluorescence on rat kidney (Kit ORG 874).

— Anti-smooth muscle antibodies are the standard diagnostic markers for autoimmune hepatitis. Anti-actin antibodies, a subgroup of smooth muscle antibodies, are found in the sera of 52–85 % of patients with active chronic autoimmune hepatitis and 22 % of patients with primary biliary cirrhosis.

SYSTEM	IFT KIT
TESTS	96 (12 slides, 8 wells each)
CODE	ORG 877
COMPONENTS	Anti-Smooth Muscle Antibodies (ASMA) slides; ASMA Positive Control; FITC Evans (R); Mounting Medium; PBS (10x); Negative Control; Blotting Paper 1x8;

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— Tissue sections of rat liver, kidney (cortex, medulla), and stomach are the substrate on the slide. They detect autoantibodies in serum or plasma against the mitochondria (AMA), smooth muscle antigens (ASMA), parietal cells of the gastric mucosa (APCA), liver-kidney microsomes (Anti-LKM) and other structures. The different types of antibodies are associated with characteristic fluorescence patterns on tissue sections.

— The presence of AMA is associated with primary biliary cirrhosis in more than 95 % of patients. ASMA are found in the sera of 52–85 % of patients with chronic autoimmune hepatitis and 22 % of patients with primary biliary cirrhosis. Up to 90 % of patients with pernicious anaemia produce APCA. Anti-LKM-1 antibodies are specific markers for autoimmune hepatitis type II.

<b>SYSTEM</b>	<b>IFT KIT</b>
<b>TESTS</b>	<b>96 (12 slides, 8 wells each)</b>
<b>CODE</b>	<b>ORG 874</b>
COMPONENTS	Autoantibodies RL/RK/RS slides; AMA Positive Control; ANA-RL Positive Control; ASMA Positive Control; FITC/Evans (R); Mounting Medium; PBS (10×); Negative Control; Blotting Paper 1×8

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— The Random Access Analyser, which is based on SMC® Technology developed by ORGENTEC Diagnostika GmbH, represents a new dimension in the field of autoimmune diagnostics and in infection serology. It provides a fully automated process for detecting antibodies with a level of flexibility never available before.

— With this system, users are now able to perform diagnostic tests in a way that is specially adapted to the individual needs of their own particular environment. Each Alegria® Test Strip contains a complete set of reagents for determining single samples. Each patient sample can be handled on an individual basis and at the lowest cost, thus liberating new resources for the benefit of the patient.

— This current development is just another milestone in the long success story of ORGENTEC Diagnostika. Thanks to its compliance with high quality requirements, the device fits seamlessly into product lines already available on the market. The innovative concept of Alegria® has got already several important design awards.

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— Software for fast and easy uni- or bi-directional connection to external laboratory information systems (LIS) via LIS standard protocols.

#### Alegria® Negative Control **ORG 311**

— 24 Test Strips for 24 individual determinations  
 — The Alegria® Negative Control is an ELISA test system intended for quality control in the Alegria® instrument. Each barcoded test strips is equipped with a complete set of reagents, including enzyme conjugate, enzyme substrate, negative sample and a test-specific control. The test result is displayed in [U/ml].

#### Alegria® Positive Control **ORG 310**

— 24 Test Strips for 24 individual determinations  
 — The Alegria® Positive Control is an ELISA test system intended for quality control in the Alegria® instrument. Each barcoded test strip is equipped with a complete set of reagents, including enzyme conjugate, enzyme substrate, positive sample and a test-specific control. The test result is displayed in [U/ml].

#### Flush Routine Solution

<b>1 x 20 ml</b>	<b>ORG 305-01</b>
<b>12 x 20 ml</b>	<b>ORG 305-12</b>

— Alegria® Flush Routine Solution is intended for use in the monthly rinse cycle of the Alegria® instrument.

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## Management of samples and reagents

- Each of the three carrier modules (SysTray A, B and C) can be loaded with ten Test Strips /patient samples.
- Up to 30 different samples with individual parameters for autoimmune or infectious disease diagnostics can be examined at the same time in one operational step.
- Serum, plasma, stool and cerebrospinal fluid samples can be analysed simultaneously in the same test run.
- Each Alegria® Test Strip has two cavities for handling samples, two incubation cavities for patient sample and a cut-off control, as well as a complete set of reagents consisting of a sample buffer, a conjugate, a substrate and an internal calibrator control.
- All reagents are packed under airtight conditions; they are not opened until immediately prior to being used, when they are directly inserted into the test in the system. This eliminates any external influence or contamination.

## SMC® Technology

- Using SMC® Technology, all test information is passed on directly to the system by means of specific barcodes.
- The SMC® coding system contains information on product names, a batch-specific standard curve, threshold values for the internal cut-off control, a recalculation algorithm feature, as well as information on production and expiry dates of individual test strips.
- Expired test strips are not accepted by the system.

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### Incubation unit

- Inside the thermally regulated incubation unit, all processes run under strictly standardised conditions.
- The system cycle is unable to start until predefined climatic conditions are reached.
- Provisions for reloading other system trays (SysTray A, B, C) can be made in advance while another test procedure is still running, thus resulting in maximum capacity use in terms of flexibility and high sample throughput.
- When system requirements are reached inside the incubation unit, test-specific parameters for each reagent strip are selected by means of SMC® Technology and assigned to the relevant sample.

### Barcode management

- The external barcode reader allocates individual patient samples to their own special location within the system. Together with the test information shown by SMC® Technology, specific evaluations are created for each patient, thus making it impossible to confuse results.

### Integrated printer

- After a complete analysis all test data corresponding to the patients tested are printed out by the integrated printer.

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### Software link-up and data management

- Overall performance and data management are both supported by Windows®, thus enabling a quick and uncomplicated link-up to external computer systems via already installed standard protocols.
- Alegria® can be programmed and adapted at any time by system technicians via Windows® based maintenance and system software.

### Touchscreen

- Thanks to its easy-to-use software, all system functions can be easily operated and controlled via a user-friendly touchscreen.

### System liquid supply

- There are four different symbol coded reagent containers located in the device (2 × wash buffers – both 1 litre; 1 × System Fluid – 2.5 litres; 1 × liquid waste container – 5 litres)
- All containers are regulated by a system of separate precision pumps, and equipped with liquid level sensors.

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### An overview of the system's most significant advantages

- Random Access System with SMC® Technology
- SMC® Technology – Sensotronic Memorized Calibration
- Robust and reproducible operating conditions, thanks to temperature controlled incubation unit
- Maximum flexibility in diagnosis of autoimmune and infectious diseases
- Over 100 test systems available
- Up to 30 different results in just 90 minutes
- Free and flexible choice of up to 30 patients, with individual results in just one operation
- Every sample is validated through its own control

## All Products A–Z

### ALEGRIA

Alegria <sup>®</sup> Instrument . . . . .	ORG 300-00	Anti-dsDNA IgG . . . . .	ORG 204G-24
Alegria <sup>®</sup> Negative Control . . . . .	ORG 311	Anti-dsDNA IgM . . . . .	ORG 204M-24
Alegria <sup>®</sup> Positive Control . . . . .	ORG 310	Anti-dsDNA Screen . . . . .	ORG 204S-24
Flush Routine Solution . . . . .	ORG 305-01	Anti-EBV (EBNA-1) IgG . . . . .	ORG 902
Flush Routine Solution . . . . .	ORG 305-12	Anti-EBV (VCA) IgG . . . . .	ORG 901G
25-OH Vitamin D <sub>3</sub> /D <sub>2</sub> . . . . .	ORG 270	Anti-EBV (VCA) IgG Liquor . . . . .	ORG 901GL
25-OH Vitamin D <sub>3</sub> /D <sub>2</sub> Control Set . . . . .	ORG 271	Anti-EBV (VCA) IgM Abs. . . . .	ORG 901MX
Alpha-1-Antitrypsin . . . . .	ORG 286	Anti-EBV (ZEBRA) IgM. . . . .	ORG 910
AMA-M2 . . . . .	ORG 216-24	Anti-Elastase . . . . .	ORG 224-24
ANA Detect. . . . .	ORG 200-24	Anti-GBM. . . . .	ORG 250-24
ANCAscreen <i>hs (high sensitive)</i> . . . . .	ORG 389	Anti-Gliadin IgA. . . . .	ORG 234A-24
Anti-alpha-Fodrin IgA . . . . .	ORG 242A-24	Anti-Gliadin IgG . . . . .	ORG 234G-24
Anti-alpha-Fodrin IgG . . . . .	ORG 242G-24	Anti-Gliadin Screen . . . . .	ORG 234S-24
Anti-Annexin V IgG . . . . .	ORG 243G-24	Anti-gp210 . . . . .	ORG 256
Anti-Annexin V IgM. . . . .	ORG 243M-24	Anti-Helicobacter pylori IgA . . . . .	ORG 917A
Anti-B. pertussis Toxin IgA. . . . .	ORG 916A	Anti-Helicobacter pylori IgG . . . . .	ORG 917G
Anti-B. pertussis Toxin IgG . . . . .	ORG 916G	Anti-Hepatitis E-Virus . . . . .	ORG 921G
Anti-beta-2-Glycoprotein I IgA . . . . .	ORG 221A-24	Anti-Hepatitis E-Virus Abs. . . . .	ORG 921MX
Anti-beta-2-Glycoprotein I IgG . . . . .	ORG 221G-24	Anti-Histone . . . . .	ORG 207-24
Anti-beta-2-Glycoprotein I IgM . . . . .	ORG 221M-24	Anti-HSV-1 IgG . . . . .	ORG 903G
Anti-beta-2-Glycoprotein I Screen . . . . .	ORG 221S-24	Anti-HSV-1 IgM Abs. . . . .	ORG 903MX
Anti-Borrelia IgG . . . . .	ORG 911G	Anti-HSV-1/2 IgG . . . . .	ORG 905G
Anti-Borrelia IgG Liquor. . . . .	ORG 911GL	Anti-HSV-1/2 IgG Liquor. . . . .	ORG 905GL
Anti-Borrelia IgM Abs. . . . .	ORG 911MX	Anti-HSV-1/2 IgM Abs. . . . .	ORG 905MX
Anti-Borrelia IgM Liquor. . . . .	ORG 911ML	Anti-HSV-2 IgG . . . . .	ORG 904G
Anti-BPI. . . . .	ORG 223-24	Anti-HSV-2 IgM Abs. . . . .	ORG 904MX
Anti-C1q . . . . .	ORG 249-24	Anti-Insulin . . . . .	ORG 220-24
Anti-Cardiolipin IgA. . . . .	ORG 215A-24	Anti-Intrinsic Factor . . . . .	ORG 247-24
Anti-Cardiolipin IgG. . . . .	ORG 215G-24	Anti-Jo-1 . . . . .	ORG 213-24
Anti-Cardiolipin IgM . . . . .	ORG 215M-24	Anti-Lactoferrin . . . . .	ORG 227-24
Anti-Cardiolipin Screen . . . . .	ORG 215S-24	Anti-LC1 . . . . .	ORG 258
Anti-Cathepsin G . . . . .	ORG 225-24	Anti-LKM-1 . . . . .	ORG 253
Anti-CCP <i>hs (high sensitive)</i> <sup>®</sup> . . . . .	ORG 301	Anti-Lysozyme . . . . .	ORG 226-24
Anti-Centromere B . . . . .	ORG 233-24	Anti-MCV <sup>®</sup> . . . . .	ORG 248-24
Anti-Chlamydia pneumoniae IgA . . . . .	ORG 907A	Anti-Measles Virus IgG . . . . .	ORG 909G
Anti-Chlamydia pneumoniae IgG . . . . .	ORG 907G	Anti-Measles Virus IgG Liquor . . . . .	ORG 909GL
Anti-Chlamydia pneumoniae IgM Abs. . . . .	ORG 907MX	Anti-Measles Virus IgM Abs. . . . .	ORG 909MX
Anti-Chlamydia trachomatis IgA. . . . .	ORG 906A	Anti-MPO. . . . .	ORG 219-24
Anti-Chlamydia trachomatis IgG. . . . .	ORG 906G	Anti-Mumps Virus IgG . . . . .	ORG 908G
Anti-Chlamydia trachomatis IgM Abs. . . . .	ORG 906MX	Anti-Mumps Virus IgM Abs. . . . .	ORG 908MX
Anti-DGP IgA . . . . .	ORG 251A	Anti-Mycoplasma pneumoniae IgA . . . . .	ORG 915A
Anti-DGP IgG . . . . .	ORG 251G	Anti-Mycoplasma pneumoniae IgG . . . . .	ORG 915G
Anti-DGP Screen . . . . .	ORG 251S	Anti-Mycoplasma pneumoniae IgM Abs.. . . .	ORG 915MX
Anti-dsDNA IgA. . . . .	ORG 204A-24	Anti-Nucleosome . . . . .	ORG 228-24
		Anti-Parietal Cell . . . . .	ORG 231-24
		Anti-Parvovirus B19 IgG. . . . .	ORG 912G
		Anti-Parvovirus B19 IgM Abs. . . . .	ORG 912MX

Anti-PR3 hs . . . . .	ORG 618
Anti-Prothrombin IgA . . . . .	ORG 541A
Anti-Prothrombin IgG/IgM . . . . .	ORG 541
Anti-Prothrombin Screen . . . . .	ORG 541S
Anti-Rib-P . . . . .	ORG 517
Anti-RNP-70 . . . . .	ORG 632
Anti-RNP/Sm . . . . .	ORG 511
Anti-Scl-70 . . . . .	ORG 512
Anti-Sm . . . . .	ORG 510
Anti-SS-A . . . . .	ORG 508
Anti-SS-A 52 . . . . .	ORG 652
Anti-SS-A 60 . . . . .	ORG 660
Anti-SS-B . . . . .	ORG 509
Anti-ssDNA . . . . .	ORG 605
Anti-TG . . . . .	ORG 502
Anti-Tissue-Transglutaminase IgA . . . . .	ORG 540A
Anti-Tissue-Transglutaminase IgG . . . . .	ORG 540G
Anti-Tissue-Transglutaminase Screen . . . . .	ORG 540S
Anti-TPO . . . . .	ORG 503
ASCA IgG/IgA . . . . .	ORG 545
beta-2-Microglobulin . . . . .	ORG 5BM
Calprotectin . . . . .	ORG 580
ENA-4-Profile . . . . .	ORG 544
ENA-6-Profile . . . . .	ORG 546
ENACombi . . . . .	ORG 514
ENAScreen . . . . .	ORG 506
Ferritin . . . . .	ORG 5FE
Micro-Albumin . . . . .	ORG 5MA
Rheumatoid Factor IgA . . . . .	ORG 522A
Rheumatoid Factor IgG . . . . .	ORG 522G
Rheumatoid Factor IgM . . . . .	ORG 522M
Rheumatoid Factor Screen . . . . .	ORG 522S
Thyroglobulin . . . . .	ORG 5TG

#### IMMUNOBLOT

ANA-9-Line . . . . .	ORG 710-08, ORG 710-16
ANCA-3-Line . . . . .	ORG 789-08, ORG 789-16
Gastro-5-Line . . . . .	ORG 740-08, ORG 740-16
Liver-9-Line 2nd Generation . . . . .	ORG 722-08, ORG 722-16
Myositis plus . . . . .	ORG 760-08, ORG 760-16
Nucleo-9-Line . . . . .	ORG 711-08, ORG 711-16

#### RAPID TEST

rheumachec <sup>®</sup> . . . . .	ORG 170-01
rheumachec <sup>®</sup> . . . . .	ORG 170-10

#### IMMUNOFLUORESCENCE

Anti-Endomysium Antibodies (AEA) – Kit . . . . .	ORG 872
Anti-nDNA Antibodies (nDNA) – Kit . . . . .	ORG 871
Anti-Neutrophil Cytoplasmic Antibodies (ANCA-Ethanol) – Kit . . . . .	ORG 891
Anti-Neutrophil Cytoplasmic Antibodies (ANCA-Formalin) – Kit . . . . .	ORG 892
Anti-Nuclear Antibodies HEp-2 (ANA-HEp-2) – Kit . . . . .	ORG 870
Anti-Smooth Muscle Antibodies (ASMA) – Kit . . . . .	ORG 877
Autoantibodies RL/RK/RS – Kit . . . . .	ORG 874

## Imprint

**ORGENTEC** Diagnostika GmbH

Carl-Zeiss-Straße 49-51

55129 Mainz, Germany

Phone +49(0)61 31/92 58-0

Fax +49(0)61 31/92 58-58

[orgentec@orgentec.com](mailto:orgentec@orgentec.com)

[www.orgentec.com](http://www.orgentec.com)

## Design

Klöppinger Identity, Darmstadt, Germany

[www.kloepfinger-identity.de](http://www.kloepfinger-identity.de)

## Art

© Roman Schmelter, Mainz, Germany

[www.galerie-schmelter.de](http://www.galerie-schmelter.de)

Supported by:

Federal Ministry  
of Economics  
and Technology

Development of many of the ORGENTEC-Tests  
for infectious disease diagnostics was supported  
by grants of the German Federal Ministry of  
Economics and Technology

**ORGENTEC** Diagnostika GmbH  
Carl-Zeiss-Straße 49-51  
55129 Mainz  
Germany

Phone +49 (0) 61 31/92 58-0  
Fax +49 (0) 61 31/92 58-58

[orgentec@orgentec.com](mailto:orgentec@orgentec.com)  
[www.orgentec.com](http://www.orgentec.com)